

Алматы (727)345-47-04
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Россия +7(495)268-04-70

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Казахстан +7(727)345-47-04

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Беларусь +(375)257-127-884

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саранск (8342)22-96-24
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Узбекистан +998(71)205-18-59

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Улан-Удэ (3012)59-97-51
Ульяновск (8422)24-23-59
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Киргизия +996(312)96-26-47

<https://becton.nt-rt.ru> || bog@nt-rt.ru

The goal is simple: the best possible outcome

BD regional anesthesia—products you can trust, from a partner you can rely on



A lot goes into an anesthesia procedure

Every patient is unique, requiring the approach that's right for them—and the right tools for the procedure. Ultimately, successful anesthesia outcomes depend upon the skill and experience of the anesthesiologist. But we're proud that so many of them trust BD products to help them manage their patients' pain to achieve better outcomes.



100 years

Of creating innovative, reliable, high-quality regional anesthesia products



A robust portfolio

Available in multiple configurations of needles, kits, sets and trays



The BD standard

Precision engineering and high quality for BD needles and syringes in every tray



Regional is recommended

Several studies **recommend regional anesthesia techniques over general anesthesia** when patient appropriate.



Superior pain relief and improvement in post-operative pain^{*1,3}



Leading to earlier, post-surgery ambulation^{*1,2,4}



Fewer in-hospital complications^{*5}



Shorter length of stay^{*6}

Spinal anesthesia

All of spinal anesthesia from one supplier—BD offers the BD® Whitacre and BD® Quincke Spinal Needles anesthesiologists want for their quality and performance, plus additional components they need to conduct the procedure. Our spinal anesthesia trays are available in both standard and custom so anesthesiologists can choose components that cater to their specific needs.



PDPH is “the most common serious complication”

Postdural puncture headache (PDPH) is the most common serious complication resulting from spinal and epidural anesthesia and lumbar punctures.⁷ PDPH is generally thought to be the result of decreased cerebrospinal fluid (CSF) pressure due to CSF leakage through the dura puncture hole.^{7,8} The size and type of spinal needles are two of the most important factors for decreasing the risk of PDPH.⁷



BD® Whitacre Pencil-Point Spinal Needle



BD® Quincke Cutting Tip Spinal Needle

Studies show that pencil-point needles result in a five-fold reduction in CSF loss and PDPH when compared to a cutting needle. The difference is significant. Similarly, using a narrow gauge 25 G needle results in a six-fold reduction in CSF loss and PDPH compared to a 22 G needle.⁷ Using pencil-point needles has also been shown to reduce post-lumbar puncture headache (PLPH) in lumbar puncture and myelogram patients.⁹


BD® Spinal Needles

BD® Whitacre and BD® Quincke Spinal Needles are selected by clinicians around the world. Their hubs are designed to enable anesthesiologists to **visualize CSF for proper needle placement**, and the **manufacturing quality, reliability and consistency** of BD products can give them **peace of mind**.




Proven experience

With the BD® Whitacre Spinal Needle supported by over 25 peer-reviewed publications



Enhanced precision

With BD® Spinal Needles' hubs designed to enable visualization of CSF, aiding proper needle placement



From a trusted partner

With a range of products and configurations to meet anesthesia needs throughout your facility

BD® Spinal Needles

BD® Whitacre Spinal Needles

Supported by over 25 peer-reviewed publications, the BD® Whitacre Spinal Needle is the pioneer of pencil-point spinal needles—designed to deliver enhanced performance, increased reliability¹⁰ and support for anesthesiologists in conducting successful spinal anesthesia procedures.

 **Distinct tactile "click" or "pop"** as the pencil point penetrates the dura, alerting the clinician of proper placement

 Associated with significantly **lower incidence of blood patch rate** post-lumbar puncture¹⁴

 Designed to provide **rapid CSF detection**^{11,12} and **reduced CSF loss**¹³

 Blunt, pencil-point tip is designed to spread tissue without cutting the tissue, potentially **lowering the risk of PDPH**¹⁵

Available in 8 sizes from 3 ½ in. to 5 in. and in gauges from 22 G to 27 G. Some BD® Whitacre Spinal Needles have thin walls, designed to maximize flow through needles of the same gauge.



Rollover points to see features.



BD® Spinal Needles

BD® Quincke Spinal Needles

BD® Quincke Spinal Needles have a unique key/slot arrangement of stylet and cannula hubs that facilitates proper needle bevel orientation.

 Lesions with BD® Quincke Spinal Needles are shown to result in a clean cut opening in the dura

 Hub provides visualization of cerebrospinal fluid (CSF)

 Wide range of needle configurations to meet varied clinical needs

Available in 22 sizes with lengths ranging from 1 in. to 7 in. and gauges ranging from 18 G to 27 G.



Rollover points to see features.




BD® Spinal Needles

BD® Spinal Anesthesia Trays

The pioneer in spinal needles, now with the quality and convenience of a BD® Spinal Anesthesia Tray

BD® Spinal Anesthesia Trays are designed for successful blocks—from preparation and placement to delivery. Available in standard, drugless and safety configurations, BD® Spinal Anesthesia Trays provide trusted needles, syringes and safe medication preparation all in one solution.




Rollover numbers to highlight in the tray.

Epidural anesthesia


A full epidural solution that helps you perform the procedure, from start to finish

The BD portfolio for epidural anesthesia—needles, syringes, infusion pumps, sets, accessories and trays—contains everything anesthesiologists need to **perform an epidural block**. Our high-quality products are designed with features that **enhance procedural control**.



Long-lasting pain management

With SPIROL® Catheters designed to reduce the risk of kinking, helping enhance patient comfort¹⁶



Facilitates clinical control

With ergonomic grips on BD Perisafe™ Weiss Epidural Needles and enhanced tactile precision on BD Epilor™ Glass Loss-of-Resistance (LOR) Syringes



From a trusted partner

With a range of products and configurations to meet anesthesia needs throughout your facility

BD® Epidural Needles

BD® Tuohy Epidural Needles

BD® Tuohy Epidural Needles feature a long, curved bevel with distal tip designed to safely enter the epidural space and facilitate catheter introduction, while a smooth proximal bevel helps reduce the risk of cutting a receding catheter.¹⁵



Both the BD® Weiss and the BD® Tuohy Epidural Needles have **calibrated needle-depth markings** to illustrate depth. So whether anesthesiologists are looking for tactile feedback or penetration characteristics, BD has a range of epidural needles to provide the control that meets their needs.



Rollover points to see features.



BD® Epidural Needles

BD® Weiss Epidural Needles

The fixed wings of the BD® Weiss Epidural Needle provide an ergonomic grip for anesthesiologists, helping them to easily guide the needle inside the spine of a patient.



Both the BD® Weiss and the BD® Tuohy Epidural Needles have **calibrated needle-depth markings** to illustrate depth. So whether anesthesiologists are looking for tactile feedback or penetration characteristics, BD has a range of epidural needles to provide the control that meets their needs.



Rollover points to see features.



BD® Epidural Needles

BD® Epidural Anesthesia Trays


Available in standard, specialty and drugless configurations for both single-shot and continuous epidural anesthesia, the BD portfolio of epidural trays contains needles, LOR syringes, accessories and everything anesthesiologists need to perform an epidural procedure.



Rollover numbers to highlight in the tray.

BD® Epidural Needles

BD® Epidural Components



BD® Epidural Catheters

The SPIROL® Catheter has the perfect balance of strength and flexibility for reliably delivering medication where it is most needed. From the soft flexible tip to the firm, kink-resistant body, the SPIROL® Catheter is designed to deliver a new measure of fault-free performance.¹⁶ BD Perisafe™ Nylon Catheters are available in both closed-end (CE) and open-end (OE) multi-port tip configurations.




BD® Loss-of-Resistance (LOR) Syringes

BD Epilor™ LOR Syringes facilitate the “loss of resistance” technique for identifying the epidural space by reducing subjectivity when locating this space and the potential for complications when administering epidural anesthesia to patients. BD also offers BD Luer-Lok™ and luer slip LOR syringes.

Combined spinal and epidural (CSE) anesthesia

It's the best of both worlds, in one convenient solution

BD offers the BD® Whitacre Spinal Needle and the BD® Weiss Epidural Needle, plus all the other components anesthesiologists need to perform CSE procedures. BD has a **robust CSE portfolio** with a variety of **matching CSE needle sets, catheters and LOR syringes**.




Enhanced experience

With SPIROL® Catheters designed to resist kinking,¹⁶ and the BD® Whitacre Spinal Needle supported by over 25 peer-reviewed publications



Greater convenience

With a matched CSE set, containing the BD® Whitacre Spinal Needle



From a trusted partner

With a range of products and configurations to meet anesthesia needs throughout your facility

Combined spinal and epidural (CSE) anesthesia

BD® CSE Trays

Available in standard, drugless, specialty and safety configurations, BD® CSE Trays deliver four benefits in one place.



Rollover numbers to highlight in the tray.

Combined spinal and epidural (CSE) anesthesia

BD® CSE Sets

BD offers kits and sets designed specifically for CSE procedures. They contain the world-class BD® Whitacre Spinal Needle with a hub designed to fit with the BD® Weiss Epidural Needle. The fixed wings of the BD® Weiss Epidural Needle provide an ergonomic grip for anesthesiologists, helping enable control during insertions.



Appendix:

Product SKU List

BD® Spinal Needles

BD® Whitacre Spinal Needles			
Order no.	Description	Order no.	Description
405010	22 G x 3 ½ in. (89 mm)	405140	25 G x 5 in. (127 mm), TW
405133	24 G x 3 ½ in. (89 mm), TW	405079	27 G x 3 ½ in. (89 mm), TW
405138	25 G x 3 ½ in. (89 mm), TW	409443	27 G x 4 ⅞ in. (119 mm), TW
409442	25 G x 4 ⅞ in. (119 mm)	405144	27 G x 5 in. (127 mm), TW

BD® Quincke Spinal Needles			
Order no.	Description	Order no.	Description
405174	18 G x 3 in. (76 mm)	405181	22 G x 3 ½ in. (89 mm)
405184	18 G x 3 ½ in. (89 mm)	405148	22 G x 5 in. (127 mm)
408360	18 G x 6 in. (152 mm)	405149	22 G x 7 in. (178 mm)
405162	20 G x 1 ½ in. (38 mm)	400106	23 G x 3 ½ in. (89 mm)
405071	20 G x 2 ½ in. (64 mm)	405073	25 G x 1 in. (25 mm) Neonatal
405172	20 G x 3 in. (76 mm)	405078	25 G x 2 in. (51 mm)
405182	20 G x 3 ½ in. (89 mm)	405170	25 G x 3 in. (76 mm)
405211	20 G x 6 in. (152 mm)	405180	25 G x 3 ½ in. (89 mm)
405161	22 G x 1 ½ in. (38 mm)	405234	25 G x 4 ⅞ in. (119 mm)
405074	22 G x 2 ½ in. (64 mm)	405164	26 G x 3 ½ in. (89 mm)
405171	22 G x 3 in. (76 mm)	405081	27 G x 3 ½ in. (89 mm)

Introducer needles	
Order no.	Description
405160	20 G x 1 ¼ in. (32 mm)

405160 - For use with 25 G and 27 G needle or smaller.

BD needle gauge color chart as per ISO 6009 standards	16 G	17 G	18 G	20 G	21 G	22 G	23 G	24 G	25 G	26 G	27 G
---	------	------	------	------	------	------	------	------	------	------	------

BD® Spinal Anesthesia Trays

		Trays	400866	400867	400868	405609	405621	405632	405663	405637	405645	405646	405652	405658	405671	405672	405673	405674	405675	405699	405707	405709	405718	405721	405720	405725	405741	405723	Drugless 405727	Safety 406704	Support 406060
Procedure		Whitacre Spinal Needles																													
Procedure needles		22 G x 3 ½ in. BD® Whitacre Spinal Needle, TW																													
		24 G x 3 ½ in. BD® Whitacre Spinal Needle, TW																													
		25 G x 3 ½ in. BD® Whitacre Spinal Needle, TW																													
		27 G x 3 ½ in. BD® Whitacre Spinal Needle, TW																													
		Quincke Spinal Needles																													
Procedure needles		22 G x 3 ½ in. BD® Quincke Spinal Needle																													
		23 G x 3 ½ in. BD® Quincke Spinal Needle																													
		25 G x 3 ½ in. BD® Quincke Spinal Needle																													
		26 G x 3 ½ in. BD® Quincke Spinal Needle																													
		Sprotte® Needles																													
Introducer needles		24 G x 3 ½ in. Sprotte® Needle																													
		18 G x 1 ½ in. introducer needle																													
		20 G x 1 ½ in. introducer needle																													
Procedure syringes		5 mL BD Glaspak® Luer Slip Tip Syringe																													
		5 mL BD® Luer Slip Tip Syringe																													
		5 mL BD Luer-Lok® Syringe																													
Drugs		Bupivacaine HCL (75%) with dextrose (8.25%), 2 mL (Marcaine®)																													
Drugs		Lidocaine HCL (1%), 5 mL																													
Medication preparation																															
Medication preparation components		18 G x 1 ½ in. blunt filter needle (5 micron filter)																													
		18 G x 1 ½ in. needle																													
		1 ½ in. filter straw (5 micron filter)																													
Hypodermic injection		Syringe label																													
		Needles																													
		21 G x 1 ½ in. needle																													
Hypodermic injection components		22 G x 1 ½ in. needle																													
		25 G x 1 in. needle																													
		25 G x 1 ½ in. needle																													
		27 G x 1 ½ in. needle																													
		22 G x 1 ½ in. BD SafetyGlide® Hypodermic Needle																													
Hypodermic injection components		25 G x 1 in. BD SafetyGlide® Hypodermic Needle																													
		Syringes																													
		3 mL BD Luer-Lok® Syringe																													
Skin preparation		3 mL BD® Luer Slip Syringe																													
		3 mL BD Luer-Lok® Syringe																													
		3 mL BD Luer-Lok® Syringe with 25 G x 1 in. needle																													
Skin preparation components		Povidone-iodine topical solution, ¼ oz.																													
		Prep sponge-sticks																													
		Drapes																													
Skin preparation components		Blue paper fenestrated drape (3 x 4 in. square, 27 x 17 ½ in.)																													
		Clear plastic fenestrated drape (4 in. round, 32 x 24 in.)																													
		Blue paper fenestrated drape (4 in. round, 30 x 26 in.)																													
Procedure accessories																															
Gauze		Gauze pads																													
Bandage		Adhesive bandage																													
Paper towel		Paper towel																													
Sharps safety		Sharps stick pad																													
		Point-Lok® Sharps Safety Device																													

BD® Epidural Single-shot Trays

BD® Epidural Single-shot Trays		Trays	406070	406072	406073	406082	406087	406091	406096	406122	406134	406157	Support 406308	400862	406069
Procedure		Weiss Epidural Needles (fixed wings)													
Procedure needles		17 G x 3 ½ in. BD Perisafe® Weiss Epidural Needle													
		18 G x 3 ½ in. BD Perisafe® Weiss Epidural Needle													
		20 G x 3 ½ in. BD Perisafe® Weiss Epidural Needle													
		Tuohy Epidural Needles													
Medication delivery syringes		18 G x 3 ½ in. BD Perisafe® Tuohy Epidural Needle													
		20 G x 3 ½ in. BD Perisafe® Tuohy Epidural Needle													
		Plastic syringes													
		1 mL BD® Luer Slip Tip Syringe													
LOR syringes		10 mL BD® Luer Slip Tip Syringe													
		20 mL BD® Luer Slip Tip Syringe													
		10 mL BD Luer-Lok® Syringe													
		20 mL BD Luer-Lok® Syringe													
		Epilor® Glass LOR Syringes													
		5 mL BD Epilor® LOR Glass Luer Slip, metal-tip syringe													
		5 mL BD Epilor® LOR Glass Luer Slip, glass-tip syringe													
		5 mL BD Epilor® LOR Glass BD Luer-Lok®, metal-tip syringe													
Epidural accessories		Epilor® Plastic LOR Syringe													
		7 mL BD Epilor® LOR Plastic Luer Slip Tip Syringe													
		BD Perisafe® Epidural Catheter Connector													
		Cap BD Luer-Lok® male													
Drugs															
Drugs		Lidocaine HCL (1%), 5 mL													
		Lidocaine HCL (1.5%) with Epi (1:200,000), 5 mL													
		Sodium chloride (0.9%), 10 mL													
Medication preparation															
Medication preparation components		Filter needles													
		18 G x 1 ½ in. needle													
		18 G x 1 ½ in. blunt filter needle (5 micron)													
		4 in. filter straw (5 micron filter)													
Hypodermic injection		Syringe label													
		Needles													
		21 G x 1 ½ in. needle													
		22 G x 1 ½ in. needle													
Hypodermic injection components		25 G x 1 in. needle													
		25 G x 1 ½ in. needle													
		27 G x 1 ¼ in. needle													
		30 G x 1 in. needle													
		Syringes													
		3 mL BD Luer-Lok® Syringe													
		3 mL luer slip syringe													
		5 mL BD Luer-Lok® Syringe													
Skin preparation															
Skin preparation components		20 mL BD Luer-Lok® Syringe													
		Povidone-iodine topical solution, ¼ oz.													
		Prep sponge-sticks													
Procedure accessories		Drapes													
		Clear plastic fenestrated drape (4 in. round, 32 x 24 in.)													
		Blue paper fenestrated drape (4 in. round, 30 x 26 in.)													
Paper towel															
Gauze															
Sharps safety															

BD® Epidural Continuous Trays

BD® Epidural Continuous Trays		Trays	400835	400840	400844	400846	406048	406123	406050	400835	400840	400847	406015	406034	406035	406047	406048	406049	406050	406151	406013
Procedure																					
Procedure needles	Weiss Epidural Needles (fixed wings)																				
	17 G x 3 ½ in. BD Perisafe® Weiss Epidural Needle		*							*		*	*	*	*					*	
	17 G x 3 ½ in. BD Perisafe® Weiss Epidural Needle (removable wings)													*							
	18 G x 1 ½ in. needle																				
	18 G x 3 ½ in. BD Perisafe® Weiss Epidural Needle			*	*	*	*	*	*		*					*	*	*	*	*	
Tuohy Epidural Needles																					
Medication delivery syringes	17 G x 3 ½ in. BD Perisafe® Tuohy Epidural Needle																				*
	Plastic syringes																				
	20 mL BD Luer-Lok® Syringe		*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*
	20 mL luer slip tip syringe												*	*	*	*	*	*	*	*	*
	Epilor® Glass LOR Syringes						*							*	*	*	*	*	*	*	*
LOR syringes	5 mL BD Epilor® LOR Glass Luer Slip, metal-tip syringe						*						*	*	*	*	*	*	*	*	*
	5 mL BD Epilor® LOR Glass Luer Slip, glass-tip syringe						*	*					*	*	*	*	*	*	*	*	*
	5 mL BD Epilor® LOR Glass BD Luer-Lok®, metal-tip syringe		*	*		*				*	*	*		*	*	*	*	*	*	*	*
	Epilor® Plastic LOR Syringe																				
	7 mL BD Epilor® LOR Plastic Luer Slip Tip Syringe				*				*										*		
Nylon catheters																					
Catheters	19 G BD Perisafe® Nylon Catheter, CE												*	*	*						*
	20 G BD Perisafe® Nylon Catheter, CE			*	*	*	*	*	*		*	*			*	*	*	*	*	*	
	Spring-reinforced catheters																				
Epidural accessories	19 G Spinal® Spring-reinforced Catheter, CE		*							*										*	
	Flat filter (0.2 micron filter)		*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*
	BD Perisafe® Epidural Catheter Connector		*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*
	Cap BD Luer-Lok® male		*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*
	Cap BD Luer-Lok® female		*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*
Catheter support pad		*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*
Drugs																					
Drugs	Lidocaine HCL (1%), 5 mL		*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*
	Lidocaine HCL (1.5%) with Epi (1:200,000), 5 mL		*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*
	Sodium chloride (0.9%), 10 mL		*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*
Medication preparation																					
Medication preparation components	Filter needles																				
	18 G x 1 ½ in. needle						*						*	*		*			*	*	*
	18 G x 1 ½ in. blunt filter needle (5 micron)		*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*
	4 in. filter straw (5 micron filter)		*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*
Syringe label		*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*
Hypodermic injection																					
Hypodermic injection components	Needles																				
	22 G x 1 ½ in. needle		*	*	*		*	*	*	*	*		*	*	*	*	*	*	*	*	*
	25 G x 1 in. needle				*		*	*	*				*	*	*	*	*	*	*	*	*
	25 G x 1 ½ in. needle		*	*		*				*	*	*		*	*	*	*	*	*	*	*
	27 G x 1 ¼ in. needle		*	*	*					*	*	*		*	*	*	*	*	*	*	*
Syringes																					
Skin preparation	3 mL BD Luer-Lok® Syringe		*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*
	3 mL luer slip syringe												*			*	*	*	*	*	*
Skin preparation components																					
Skin preparation components	Povidone-iodine topical solution, ¼ oz.		*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*
	Prep sponge-sticks		*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*
	Drapes																				
	Clear plastic fenestrated drape (4 in. round, 32 x 24 in.)		*	*	*	*		*	*	*	*	*	*				*	*	*	*	*
Blue paper fenestrated drape (4 in. round, 30 x 26 in.)				*										*	*	*	*			*	*
Paper non-fenestrated drape		*			*					*											*
Procedure accessories																					
Paper towel	Paper towel		*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*
Gauze	Gauze pads		*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*
Sharps safety	Sharps stick pad		*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*	*

BD® Epidural Needles and Components

BD® Weiss Epidural Needles			
Order no.	Description	Order no.	Description
405194	16 G x 3 ½ in. (90 mm)	405190	18 G x 5 in. (127 mm)
405191	17 G x 3 ½ in. (90 mm)	405185	20 G x 2 in. (51 mm)
405193	17 G x 5 in. (127 mm)	405186	20 G x 3 ½ in. (90 mm)
405188	18 G x 3 ½ in. (90 mm)		

BD® Tuohy Epidural Needles			
Order no.	Description	Order no.	Description
405092	17 G x 3 ½ in. (90 mm)	405028	20 G x 3 ½ in. (90 mm)
405018	18 G x 3 ½ in. (90 mm)		

BD needle gauge color chart

16 G

17 G

18 G

20 G

21 G

22 G

23 G

24 G

25 G

26 G

27 G

BD Epilor® LOR Syringes	
Order no.	Description
Only available in trays	5 mL BD Epilor™ LOR Glass BD Luer-Lok™, metal-tip syringe
Only available in trays	5 mL BD Epilor™ LOR Glass Luer Slip Tip, glass-tip syringe
Only available in trays	10 mL BD Epilor™ LOR Glass Luer Slip, glass-tip syringe
405291	7 mL Luer-Lok™, plastic
405292	7 mL luer-slip, plastic

BD Perisafe® Nylon Catheters with Epidural Catheter Connector	
Order no.	Description
Only available in trays	19 G x 36 in. (914 mm) CE
400507	19 G x 36 in. (914 mm) OE
Only available in trays	20 G x 36 in. (914 mm) CE

BD® CSE Trays

		Trays	405811	405812	Safety 405814	Drugless 405823	405827	405828	400832
Procedure									
Spinal needles	25 G x 4 1/4 in. BD® Whitacre Spinal Needle		*	*	*	*	*	*	*
	27 G x 4 1/4 in. BD® Whitacre Spinal Needle								
Epidural needles	17 G x 3 1/2 in. BD Persafe® Weiss Epidural Needle (fixed wings)		*	*	*	*	*	*	*
	18 G x 3 1/2 in. BD Persafe® Weiss Epidural Needle (fixed wings)								
Medication delivery syringes	Plastic syringes								
	1 mL luer slip tip syringe			*					*
	5 mL luer slip tip syringe					*			
	10 mL luer slip tip syringe				*				
	20 mL BD Luer-Lok® Syringe		*	*		*	*	*	*
	Glass syringes								
LOR syringes	5 mL BD Glaspak® Luer Slip Tip Syringe		*		*		*		*
	Glass LOR syringes								
	5 mL BD Epilor® LOR Glass BD Luer-Lok®, metal-tip syringe		*	*			*	*	
	5 mL BD Epilor® LOR Glass Luer Slip, glass-tip syringe								*
	10 mL BD Epilor® LOR Glass Luer Slip, glass-tip syringe				*				
	Plastic LOR syringes								
Catheters	7 mL BD Epilor® LOR Plastic Luer Slip Tip Syringe				*				
	7 mL BD Epilor® LOR Plastic BD Luer-Lok® Syringe					*			
	Nylon catheters								
	20 G BD Persafe® Nylon Catheter, CE		*	*	*	*			
	Spring reinforced catheters								
	19 G Spinal® Spring-reinforced Catheter, CE						*	*	*
Epidural accessories	Flat filter (0.2 micron filter)		*	*			*	*	*
	BD Persafe® Epidural Catheter Connector		*	*	*	*			*
	Cap BD® Luer-Lok® female with luer cap		*	*	*	*	*	*	*
	Cap BD® Luer-Lok® male with luer cap		*	*	*	*	*	*	*
	Catheter support pad			*	*	*	*	*	*
Drugs									
Drugs	Lidocaine HCL (1%), 5 mL		*	*	*		*	*	*
	Lidocaine HCL (1.5%) with Epi (1200,000), 5 mL		*	*			*	*	*
	Sodium chloride (0.9%), 10 mL		*	*	*		*	*	*
Medication preparation									
Medication preparation components	18 G x 1 1/2 in. blunt filter needle (5 micron filter)		*	*	*	*	*	*	*
	18 G x 1 1/2 in. needle, short bevel				*	*	*	*	
	1 1/4 in. filter straw (5 micron filter)								*
	Syringe label		*	*	*	*	*	*	*
Skin preparation									
Skin preparation	Povidone-iodine topical solution, 1% oz.		*	*	*		*	*	*
	Prep sponge sticks		*	*	*	*	*	*	*
	Clear plastic fenestrated drape (4 in. round, 32 x 24 in.)		*	*	*	*	*	*	*
Hypodermic injection									
Hypodermic injection components	Needles								
	22 G x 1 1/2 in. needle		*	*		*	*		*
	25 G x 1 in. needle		*			*	*		
	25 G x 1 1/2 in. needle					*		*	*
	27 G x 1 1/2 in. needle			*					*
	22 G x 1 1/2 in. BD SafetyGlide® Hypodermic Needle				*				
	25 G x 1 in. BD SafetyGlide® Hypodermic Needle				*				
	Syringes								
Procedure accessories	3 mL BD Luer-Lok® Syringe		*	*	*	*	*	*	*
	3 mL luer slip tip syringe			*				*	
Gauze			*	*	*	*	*	*	*
Paper towel			*	*	*	*	*	*	*
Sharps safety	Sharps stick pad		*	*	*	*	*	*	*
	Point-Lok® Sharps Safety Device				*				

BD® CSE Sets

		Sets	405139	405452	405453
Procedure					
Spinal needles	BD® Whitacre Spinal Needles				
	25 G x 4 1/4 in. BD® Whitacre Spinal Needle		*		
	27 G x 4 1/4 in. BD® Whitacre Spinal Needle			*	*
Epidural needles	BD® Weiss Epidural Needles (fixed wings)				
	17 G x 3 1/2 in. BD Persafe™ Weiss Epidural Needle		*		*
	18 G x 3 1/2 in. BD Persafe™ Weiss Epidural Needle			*	

Алматы (727)345-47-04
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Россия +7(495)268-04-70

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Казахстан +7(727)345-47-04

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Беларусь +(375)257-127-884

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саранск (8342)22-96-24
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Узбекистан +998(71)205-18-59

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Улан-Удэ (3012)59-97-51
Ульяновск (8422)24-23-59
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Киргизия +996(312)96-26-47