

Наборы для внутривенного вливания

Техническое описание

По вопросам продаж и поддержки обращайтесь:

Алматы (727)345-47-04
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70


Казахстан +(727)345-47-04

Беларусь +(375)257-127-884

Узбекистан +998(71)205-18-59

Киргизия +996(312)96-26-47

эл.почта: bog@nt-rt.ru || сайт: <https://becton.nt-rt.ru/>




**BD MaxZero™ Extension Set,
minibore tubing, 1 bonded BD
MaxZero™ Needle-free
Connector, needle-free access
T-connector, 1 white slide
clamp**

SKU: MZXT9001




**BD MaxZero™ Extension Set,
pressure rated, standard bore,
1 removable BD MaxZero™
Needle-free Connector, 1 NAC
T-Connector, 1 slide clamp,
spin male luer lock**

SKU: MZXT5306




**BD MaxZero™ Extension Set,
pressure rated, standard bore,
1 bonded BD MaxZero™
Needle-free Connector, 1 NAC
T-Connector, 1 slide clamp,
spin male luer lock**

SKU: MZXT5304




**BD MaxZero™ Extension Set,
pressure rated, standard bore,
1 removable BD MaxZero™
Needle-free Connector, 1
NAC-Y, 1 white slide clamp,
spin male luer lock**

SKU: MZX5306




**BD MaxZero™ Extension Set,
pressure rated, standard bore,
1 bonded BD MaxZero™
Needle-free Connector, 1
NAC- Y, 1 white slide clamp,
spin male luer lock**

SKU: MZX5305




**MaxZero™ extension set,
minibore 1 MaxZero™
needle-free connector(s) 1.2
micron small volume filter
(0.25 sq. in low protein
binding) slide clamp(s)**

SKU: MZ9329




**MaxZero™ Extension Set,
Microbore Trifuse: Fuse 1: 1
Bonded MaxZero™ Needle-
Free Connector + White Slide
Clamp; Fuse 2 And 3: Female
Luers, White Slide Clamps**

SKU: MZ9327




**MaxZero™ Extension Set,
Minibore, Trifuse: Fuse 1: 1
Bonded MaxZero™ Needle-
Free Connector + Slide
Clamp; Fuse 2 And 3: Female
Luers, 1 Red and 1 Blue Slide
Clamp**

SKU: MZ9326




**MaxZero™ Extension Set,
Minibore, Trifuse, 1
MaxZero™ Needle-Free
Connector, 0.2 Micron Small
Volume Filter (0.25 Sq. In
Low Protein Binding), 3
Colored Slide Clamps, 2
Female Luers, Spin Male Luer
Lock, 1.2 Micron Filter**

SKU: MZ9325




**MaxZero™ Extension Set,
Minibore, Quadfuse, 2
MaxZero™ Needle-Free
Connectors, 0.2 Micron Small
Volume Filter (0.25 Sq. In
Low Protein Binding), 4 Slide
Clamps, 2 Female Luers Spin
Male Luer Lock, 1.2 Micron
Filter**

SKU: MZ9324



**MaxZero™ Extension Set,
Minibore, 2 Check Valves
Quadfuse 1 MaxZero™
Needle-Free Connector, 1.2
Micron Small Volume Filter
(0.25 Sq. In Low Protein
Binding), 4 Slide Clamps 3**





**MaxZero™ extension set,
microbore bifuse, 2 removable
MaxZero™ needle-free
connectors, 2 slide clamps, 0.2
micron filter on one lead,**

SKU: MZ9294


Female Luers Spin Male Luer Lock

SKU: MZ9322




MaxZero™ Extension Set, Microbore, Bifuse, 2 Removable MaxZero™ Needle-Free Connectors, 2 Slide Clamps, Check Valve On Only 1 Lead, 0.2 Micron Filter On Other Lead, Spin Male Luer Lock

SKU: MZ9293



MaxZero™ Extension Set, minibore quadfuse, 4 bonded MaxZero™ Needle-free Connectors, 4 check valves, 4 white slide clamps

SKU: MZ9283




MaxZero™ Extension Set Minibore, Trifuse, 3 Bonded MaxZero™ Needle-Free Connectors, 3 Check Valves, 3 Slide Clamps (Red, White, Blue), Trifurcated Male Luer Lock

SKU: MZ9281


MaxZero™ Extension Set, Minibore Bifuse, 2 Bonded MaxZero™ Needle-Free Connectors, 2 Check Valves, 2 White Slide Clamps, Spin Male Luer Lock

SKU: MZ9279




**MaxZero™ Extension Set,
Microbore Trifuse, 3 Bonded
MaxZero™ Needle-Free
Connectors, 1.2 Micron Filter,
3 Slide Clamps, Spin Male
Luer Lock**

SKU: MZ9276




**MaxZero™ Extension Set,
microbore quadfuse, 4 bonded
MaxZero™ Needle-free
Connectors, 4 slide clamps,
spin male luer lock**

SKU: MZ9275




**MaxZero™ Extension Set,
microbore trifuse, 3 bonded
MaxZero™ Needle-free
Connectors, two 0.2 micron
filters, 3 slide clamps (red,
white, blue)**

SKU: MZ9274




**MaxZero™ Extension Set,
microbore trifuse, 3 bonded
MaxZero™ Needle-free
Connectors, 0.2 micron filter,
3 slide clamps (red, white,
blue)**

SKU: MZ9273




**MaxZero™ Extension Set,
Microbore, Trifuse, 3 Bonded
MaxZero™ Needle-Free
Connectors, 0.2 Micron Filter,
3 Slide Clamps, Spin Male
Luer Lock**

SKU: MZ9272



**MaxZero™ Extension Set,
Microbore, Bifuse, 2 Bonded
MaxZero™ Needle-Free
Connectors, 0.2 Micron Filter,
2 Slide Clamps, Spin Male
Luer Lock**

SKU: MZ9271



**MaxZero™ Extension Set, 3
bonded MaxZero™ Needle-
free Connectors, 0.2 micron
filter, 3 slide clamps (red,
white, blue), spin male luer
lock**


SKU: MZ9270



**BD MaxZero™ Extension Set
PVC microbore, 1 BD
MaxZero™ Needle-free
Connector(s), bonded 1 slide
clamp(s) RLL, length 152 cm
60 in.**


64-IN MICROBORE EXT. SET W/
MAXPLUS

SKU: MZ9267




**MaxZero™ Extension Set,
microbore trifuse, 3 bonded
MaxZero™ Needle-free
Connectors, 3 slide clamps
(red, white, blue), spin male
luer lock**

SKU: MZ9266




**MaxZero™ Extension Set,
microbore bifuse, 2 bonded
MaxZero™ Needle-free
Connectors, 2 slide clamps,
spin male luer lock**

SKU: MZ9265




**MaxZero™ Extension Set,
microbore, 1 bonded
MaxZero™ Needle-free
Connector, 0.2 micron filter, 1
slide clamp, spin male luer
lock**

SKU: MZ9226




**BD MaxZero™ Extension Set,
pressure rated, standard bore,
1 removable BD MaxZero™
Needle-free Connector, slide
clamp, spin male luer lock**

SKU: MZ5310



**BD MaxZero™ Extension Set,
pressure rated, standard bore,
1 bonded BD MaxZero™
Needle-free Connector, slide
clamp, spin male luer lock**

SKU: MZ5309



**BD MaxZero™ Extension Set,
pressure rated, minibore,
bifuse, 2 bonded BD
MaxZero™ Needle-free
Connectors, 2 slide clamps,
spin male luer lock**

SKU: MZ5307



**BD MaxZero™ Extension Set,
pressure rated, minibore, 1
bonded BD MaxZero™
Needle-free Connector, no
clamp, spin male luer lock**

SKU: MZ5306



**BD MaxZero™ Extension Set,
pressure rated, standard bore,
1 bonded BDMaxZero™
Needle-free Connector, no
clamp, spin male luer lock**

SKU: MZ5305




**BD MaxZero™ Extension Set,
pressure rated, standard bore,
1 removable BD MaxZero™
Needle-free Connector, pinch
clamp, spin male luer lock**

SKU: MZ5304




**BD MaxZero™ Extension Set,
pressure rated, standard bore,
1 bonded BD MaxZero™
Needle-free Connector, pinch
clamp, spin male luer lock**

SKU: MZ5303




**MaxZero™ Extension Set,
microbore trifuse, 3 bonded
MaxZero™ Needle-free
Connectors, 3 check valves, 3
white slide clamps**

SKU: MZ9278




**MaxZero™ Extension Set,
microbore, 1 bonded
MaxZero™ Needle-free
Connector, slide clamp, spin
male luer lock**


SKU: MZ9277




**BD MaxZero™ Extension Set,
pressure rated, minibore, 1
removable BD MaxZero™
Needle-free Connector, slide
clamp(s), spin male luer lock**
SKU: MZ5302



**BD MaxZero™ Extension Set,
pressure rated, minibore, 1
BD MaxZero™ Needle-free
Connector, slide clamp(s) spin
male luer lock**
SKU: MZ5301




**Extension set pressure rated
minibore bifuse 3 slide
clamp(s) 2 female luers
needle-free access t-connector
fixed male luer lock. Not made
with DEHP. L: 5.5 in L: 14 cm
PV: 0.53 mL content sterile
*labeled for use with low
Pressure power injectors up to
325 PSI and maximum flow
rate of 10 mL/second**



**Extension set minibore slide
clamp(s) needle-free access t-
connector male luer slip
adapter. Not made with
DEHP. L: 6 in L: 15 cm PV:
0.8 mL content sterile**
SKU: MXT1004

SKU: MXT1006



Extension set pressure rated minibore purple striped slide clamp(s) needle-free access t-connector spin male luer lock. Not made with DEHP. L: 6 in L: 15 cm PV: 0.8 mL content sterile *labeled for use with low Pressure power injectors up to 325 PSI and maximum flow rate of 10 mL/second

SKU: MXT1003




Extension set check valve 2 NAC-y needle-free valve(s) 8", 41" from spin male luer lock. Not made with DEHP. L: 50 in L: 127 cm PV: 8.4 mL content sterile


EXT SET STD BORE 50" 2NFC-MP

BCV NO CLMP SPIN LL

SKU: MX9178




Extension set slide clamp(s) 1 NAC-y needle-free valve(s) 6" from spin male luer lock. Not made with DEHP. L: 32 in L:



Extension set 2 pinch clamp(s) 2 male luers 2 NAC-y needle-free valve(s) 8", 31" from spin male luer lock. Not made with


**81 cm PV: 3.5 mL content
sterile**

SKU: MX9175



**Extension set 1.2 micron filter
(1.6 sq. in low protein
binding) pinch clamp(s) 1
NAC-y needle-free valve(s) 6"
from spin male luer lock. Not
made with DEHP. L: 16 in L:
41 cm PV: 5.8 mL content
sterile**


SKU: MX9171



**Extension set 2 pinch clamp(s)
3 NAC-y needle-free valve(s)
7.5", 10.5", 13.5" from spin
male luer lock. Not made with**


**DEHP. L: 40 in L: 102 cm PV:
7.5 mL content sterile**

SKU: MX9172



**Extension set 0.2 micron filter
(1.6 sq. in low protein
binding) pinch clamp(s) 1
NAC-y needle-free valve(s) 6"
from spin male luer lock. Not
made with DEHP. L: 16 in L:
40.6 cm PV: 5.8 mL content
sterile**

SKU: MX9166



**Extension set slide clamp(s) 1
NAC-y needle-free valve(s)
1.75" from spin male luer lock
kink resistant. Not made with**


**DEHP. L: 19 in L: 48 cm PV:
3.8 mL content sterile**

SKU: MX9076



**Extension set slide clamp(s) 2
NAC-y needle-free valve(s) 7",
28" from spin male luer lock.
Not made with DEHP. L: 36 in
L: 91 cm PV: 6 mL content
sterile**

EXT SET STD BORE 40" 2NFC-MP
2PINCH CLMPS SPIN LL
SKU: MX9059



**Extension set pressure rated
minibore slide clamp(s) 1
NAC-y needle-free valve(s) 3"
from spin male luer lock. Not
made with DEHP. L: 6 in L:
15 cm PV: 0.7 mL content**


**DEHP. L: 8 in L: 20.3 cm PV:
1.2 mL content sterile**

SKU: MX9060



**Extension set slide clamp(s) 2
NAC-y needle-free valve(s)
2.5", 4.25" from fixed male
luer lock. Not made with
DEHP. L: 8.5 in L: 22 cm PV:
2 mL content sterile**

SKU: MX9037



**Extension set pressure rated
pinch clamp(s) 1 NAC-y
needle-free valve(s) 4" from
spin male luer lock. Not made
with DEHP. L: 8 in L: 20 cm
PV: 0.8 mL content sterile**

**sterile *labeled for use with
low Pressure power injectors
up to 325 PSI and maximum
flow rate of 10 mL/second**

SKU: MX5302



**Extension set removable slide
clamp(s) 2 ganged 4-way
stopcocks (off indicator) 1
NAC-y needle-free valve(s) 7"
from fixed male luer lock. Not
made with DEHP. L: 64 in L:
163 cm PV: 6.5 mL content
sterile**

EXT SET STD BORE 64" 1NFC-MP 2-
4/W STPCK REMVBLE SLIDE CLMP
FIX LL

SKU: MX4450

***labeled for use with low
Pressure power injectors up to
325 PSI and maximum flow
rate of 10 mL/second**


SKU: MX5301



**Extension set slide clamp(s) 4-
way stopcock (off indicator) 2
NAC-y needle-free valve(s) 6",
23.5" from fixed male luer
lock. Not made with DEHP. L:
40 in L: 102 cm PV: 7 mL
content sterile**

EXT SET STD BORE 40" 2NFC-MP 4/W
STPCK SLIDE CLMP FIX LL

SKU: MX4439




**MaxPlus™ extension set,
minibore 1 removable
MaxPlus™ clear needle-free
connector(s) slide clamp(s)
needle-free access t-connector**

SKU: MPXT1005-C




**MaxPlus™ extension set,
minibore, 1 removable
MaxPlus™ clear needle-free
connector(s), slide clamp(s),
needle-free access t-connector**

SKU: MPXT1004-C




**MaxPlus™ Extension Set,
pressure rated, minibore, 1
removable MaxPlus™ Clear
Needle-free Connector(s),
slide clamp(s), needle-free
access, t-connector**

SKU: MPXT1003-C




**BD MaxPlus™ Extension Set,
pressure rated, 1 removable
BD MaxPlus™ Clear Needle-
free Connector, pinch
clamp(s), needle-free access,
T-Connector**

SKU: MPXT1001-C




**MaxPlus™ Extension Set,
Check Valve, 1 MaxPlus™
Clear Needle-Free Connector,
1 NAC-y Needle-Free Valve 6"
from Spin Male Luer Lock**

SKU: MPX9108-C




**MaxPlus™ extension set, 1
MaxPlus™ clear needle-free
connector(s), 1 NAC-y needle-
free valve(s), 5.5" from spin
male luer lock, kink resistant**

SKU: MPX9102-C




**MaxPlus™ Extension Set,
Minibore, 1 MaxPlus™ Clear
Needle-Free Connector, Slide
Clamps, 1 NAC-y Needle-Free
Valve 2.25" from Spin Male
Luer Lock**

SKU: MPX9101-C




**MaxPlus™ Extension Set,
Minibore, 1 MaxPlus™
Needle-Free Connector, Slide
Clamps, 1 NAC-y Needle-Free
Valve 2" from Spin Male Luer
Lock**

SKU: MPX9101




**BD MaxPlus™ Extension Set,
pressure rated, purple striped,
1 BD MaxPlus™ Clear
Needle-free Connector, slide
clamps, 1 NAC-y Needle-free
Valve 3" from spin male luer
lock**

SKU: MPX5305-C




**MaxPlus™ extension set,
pressure rated, purple striped,
1 MaxPlus™ clear needle-free
connector(s), slide clamp(s)**

SKU: MPX5304-C




**BD MaxPlus™ Extension Set,
pressure rated, purple striped,
1 removable BD MaxPlus™
Clear Needle-free Connector,
pinch clamp(s)**

SKU: MPX5302-C




**BD MaxPlus™ Extension Set,
pressure rated, purple striped,
1 removable BD MaxPlus™
Clear Needle-free Connector,
slide clamp(s)**

SKU: MPX5301-C




**BD MaxPlus™ Extension Set,
pressure rated, purple striped,
1 BD MaxPlus™ Clear
Needle-free Connector, pinch
clamp(s)**

SKU: MPX5300-C



**MaxPlus™ extension set,
minibore, bifuse, 2 removable
MaxPlus™ clear needle-free
connector(s), 2 removable
slide clamp(s), spin male luer
lock**

SKU: MP9240-C




**MaxPlus™ Extension Set,
Bifuse, 2 Removable
MaxPlus™ Clear Needle-Free
Connectors, 2 Slide Clamps,
Spin Male Luer Lock**

SKU: MP9238-C




**MaxPlus™ Extension Set,
minibore, 1 MaxPlus™ Clear
Needle-free Connector(s),
removable slide clamp(s), spin
male luer lock**

SKU: MP9236-C



**MaxPlus™ Extension Set,
Minibore, 1 MaxPlus™
Needle-Free Connector,
Removable Slide Clamp, Spin
Male Luer Lock**

SKU: MP9229



**MaxPlus™ extension set,
minibore, 1 removable
MaxPlus™ clear needle-free
connector(s), removable slide
clamp(s) standard injection,
T-connector**


SKU: MP9221-C



Picture is representative of product type


**MaxPlus™ extension set,
trifuse, 3 removable
MaxPlus™ Clear needle-free
connector(s), 3 removable
slide clamp(s), spin male luer
lock**

SKU: MP9220-C



**MaxPlus™ Extension Set,
Minibore, Bifuse, 1
Removable MaxPlus™ Clear
Needle-Free Connector, 2
Removable Slide Clamps, Spin
Male Luer Lock**

SKU: MP9218-C




**MaxPlus™ extension set,
minibore, 1 removable
MaxPlus™ Clear needle-free
connector(s), removable slide
clamp(s), standard injection t-
connector**

SKU: MP9216-C



**MaxPlus™ extension set, 1
MaxPlus™ Clear needle-free
connector(s), pinch clamp(s),
spin male luer lock kink
resistant**

SKU: MP9213-C




**MaxPlus™ Extension Set,
Minibore, Bifuse, 2
Removable MaxPlus™
Needle-Free Connectors, 2
Removable Slide Clamps, Spin
Male Luer Lock**

SKU: MP9027




**MaxPlus™ Extension Set,
minibore, 1 removable
MaxPlus™ Clear Needle-free
Connector(s), removable slide
clamp(s), spin male luer lock**

SKU: MP9014-C



**MaxPlus™ Extension Set,
minibore, 1 removable
MaxPlus™ Needle-free
Connector(s), removable slide
clamp(s), spin male luer lock**
SKU: MP9014




**MaxPlus™ Extension Set,
Minibore, 1 Removable
MaxPlus™ Needle-Free
Connector, Removable Slide
Clamps, Spin Male Luer Lock**
SKU: MP9004



**MaxPlus™ Extension Set, 1
MaxPlus™ Clear Needle-free
Connector(s), removable slide
clamp(s), spin male luer lock
kink resistant**
SKU: MP9002-C




**MaxPlus™ extension set, 1
MaxPlus™ needle-free
connector(s), removable slide
clamp(s), spin male luer lock
kink resistant**
SKU: MP9002



**MaxPlus™ Extension Set,
Minibore, 1 MaxPlus™
Needle-Free Connector,
Removable Slide Clamps, Spin
Male Luer Lock**

SKU: MP9001



**BD MaxPlus™ Extension Set,
pressure rated minibore,
purple striped, bifuse, 2
removable BD MaxPlus™
Clear Needle-free Connectors,
2 slide clamps, spin male luer
lock**

SKU: MP5315-C



**BD MaxPlus™ Extension Set,
pressure rated, purple striped,
1 removable BD MaxPlus™
Clear Needle-free Connector,
slide clamp(s), spin male luer
lock**

SKU: MP5314-C



**BD MaxPlus™ Extension Set,
pressure rated minibore,
purple striped 1 removable
BD MaxPlus™ Clear Needle-
free Connector, slide clamp(s)
spin male luer lock**

SKU: MP5313-C



**BD MaxPlus™ Extension Set,
pressure rated, purple striped,
1 BD MaxPlus™ Clear
Needle-free Connector, spin
male luer lock**

SKU: MP5312-C



**BD MaxPlus™ Extension Set,
pressure rated, minibore,
purple striped, 1 BD
MaxPlus™ Clear Needle-free
Connector, spin male luer lock**


SKU: MP5310-C



Picture is representative of product type

**MaxPlus™ Extension Set,
minibore, trifuse, 3
MaxPlus™ Clear Needle-free
Connector(s), 3 slide clamp(s),
spin male luer lock**

SKU: MP9233-C




**MaxPlus™ extension set,
minibore, bifuse, 2 MaxPlus™
clear needle-free connector(s),
2 removable slide clamp(s),
fixed male luer lock**

SKU: MP9232-C




**BD MaxPlus™ Extension Set,
pressure rated, purple
striped, 1 removable BD
MaxPlus™ Clear Needle-free
Connector, spin male luer lock**
SKU: MP5305-C



**BD MaxPlus™ Extension Set,
pressure rated, minibore,
purple striped, 1 BD
MaxPlus™ Clear Needle-free
Connector, slide clamp(s), spin
male luer lock**
SKU: MP5304-C



**BD MaxPlus™ Extension Set,
pressure rated, purple striped,
1 BD MaxPlus™ Clear
Needle-free Connector, pinch
clamp(s), spin male luer lock**
SKU: MP5303-C



**BD MaxPlus™ Extension Set,
pressure rated, minibore,
purple striped, 1 BD
MaxPlus™ Clear Needle-free
Connector, slide clamp(s), spin
male luer lock**
SKU: MP5302-C



**MaxPlus™ Extension Set,
Pressure Rated, Minibore,
Purple Striped, 1 MaxPlus™
Needle Free Connector, Slide
Clamp(s), Spin Male Luer
Lock**

SKU: MP5302



**BD MaxPlus™ Extension Set,
pressure rated, purple striped,
1 BD MaxPlus™ Clear
Needle-free Connector, pinch
clamp(s), spin male luer lock**

SKU: MP5301-C



**MaxPlus™ extension set,
pressure rated, purple striped,
1 MaxPlus™ needle-free
connector(s) pinch clamp(s)
spin male luer lock.**

SKU: MP5301




**MaxPlus™ Extension Set,
Minibore, 1 MaxPlus™
Needle-Free Connector,
Removable Slide Clamp(s),
Fixed Male Luer Lock**

SKU: MP2002



**Extension set 3 port manifold
(3 one way valves) spin male
luer lock. Not made with
DEHP. L: 44 in L: 112 cm PV:
7 mL content sterile**

44 IN STANDARD BORE EXTENSION
SET
SKU: MM4021



**Extension set flow regulator
slide clamp(s) 1 NAC-y
needle-free valve(s) 7" from
spin male luer lock. Not made
with DEHP. L: 14 in L: 36cm
PV: 2.6mL content sterile**

SKU: MFS141




**Extension set flow regulator
slide clamp(s) 1 NAC-y
needle-free valve(s) 7" from
spin male luer lock. Not made
with DEHP. L: 23 in L: 58.4
cm content sterile**

20 IN STANDARD BORE SET
SKU: MFS106




**Extension set pressure rated
purple striped spin male luer
lock. Not made with DEHP. L:
7.5 in L: 19 cm PV: 0.4 mL
content sterile *labeled for use
with low Pressure power
injectors up to 325 PSI and
maximum flow rate of 10
mL/second**

SKU: ME5304



**Extension set pressure rated
purple striped slide clamp(s)
spin male luer lock. Not made
with DEHP. L: 7.5 in L: 19 cm
PV: 0.4 mL content sterile
*labeled for use with low
Pressure power injectors up to
325 PSI and maximum flow
rate of 10 mL/second**

SKU: ME5303



**Extension set pressure rated
minibore purple striped slide
clamp(s) spin male luer lock.
Not made with DEHP. L: 6.5
in L: 16 cm PV: 0.2 mL
content sterile *labeled for use
with low Pressure power
injectors up to 325 PSI and
maximum flow rate of 10
mL/second**

SKU: ME5302



**Extension set pressure rated
purple striped pinch clamp(s)
spin male luer lock. Not made
with DEHP. L: 7.5 in L: 19 cm
PV: 0.4 mL content sterile
*labeled for use with low
Pressure power injectors up to
325 PSI and maximum flow
rate of 10 mL/second**

SKU: ME5301



**Extension set fixed male luer
lock. Not made with DEHP. L:
60 in L: 152 cm PV: 9.4 mL
content sterile**

EXT SET STD BORE 60" PINCH CLMP
2-PIECE LL
SKU: ME2068



Extension set microbore slide clamp(s) spin male luer lock. Not made with DEHP. L: 60 in L: 152 cm PV: 0.3 mL content sterile

EXT SET MICROBORE 62" SLIDE CLMP SPIN LL
SKU: ME2020




Extension set ultra-minibore pinch clamp(s) fixed male luer lock. Not made with DEHP. L: 60 in L: 152.4 cm PV: 1 mL content sterile

SKU: ME2017




Extension set minibore removable slide clamp(s) spin male luer lock. Not made with DEHP. L: 60 in L: 152.4 cm PV: 1.8 mL content sterile

EXT SET MIN BORE 60" REMVBLE SLIDE CLMP SPIN LL
SKU: ME2010




Extension set minibore pinch clamp(s) standard injection t-connector spin male luer lock. Not made with DEHP. L: 7 in L: 18 cm PV: 0.3 mL content sterile

SKU: ME1257




Extension set microbore trifuse 0.2 micron small volume filter (0.25 sq. in low protein binding) 4 pinch clamp(s) 3 female luer lock spin male luer lock. Not made with DEHP. L: 10 in L: 25 cm PV: 0.6 mL content sterile

SKU: ME1254



Extension set 0.2 micron filter (0.98 sq. in low protein binding) slide clamp(s) spin male luer lock. Not made with DEHP. L: 6 in L: 15 cm PV: 3.2 mL content sterile

SKU: ME1250




Extension set bifuse 2 slide clamp(s) 2 female luer lock fixed male luer lock. Not made with DEHP. L: 7 in L: 18 cm PV: 1.6 mL content sterile

SKU: ME1247




Extension set roller clamp(s) 2 male luer lock fixed male luer lock. Not made with DEHP. L: 40 in L: 102 cm PV: 6 mL content sterile

EXT SET STD BORE 40" PINCH CLMP
2MALE SPIN LL
SKU: ME1242



**Extension set minibore
removable slide clamp(s)
standard injection t-connector
male luer slip adapter. Not
made with DEHP. L: 5 in L:
13 cm PV: 0.3 mL content
sterile**

SKU: ME1233



**Extension set minibore
removable slide clamp(s)
standard injection t-connector
spin male luer lock. Not made
with DEHP. L: 5 in L: 13 cm
PV: 0.3 mL content sterile**

SKU: ME1232




**Extension set pinch clamp(s) 2
male luers fixed male luer
lock. Not made with DEHP. L:
40 in L: 102 cm PV: 6 mL
content sterile**

SKU: ME1230



**Extension set slide clamp(s)
spin male luer lock. Not made
with DEHP. L: 5 in L: 13 cm
PV: 0.7 mL content sterile**

SKU: ME1229



Extension set minibore slide clamp(s) spin male luer lock. Not made with DEHP. L: 6 in L: 15 cm PV: 0.3 mL content sterile

SKU: ME1227




Extension set minibore pinch clamp(s) spin male luer lock. Not made with DEHP. L: 30 in L: 76 cm PV: 0.4 mL content sterile

SKU: ME1226




Extension set 1.2 micron filter (1.6 sq. in low protein binding) slide clamp(s) spin male luer lock. Not made with DEHP. L: 15 in L: 38 cm PV: 4.8 mL content sterile

SKU: ME1225



Extension set removable slide clamp(s) spin male luer lock. Not made with DEHP. L: 7 in L: 18 cm PV: 1 mL content sterile

SKU: ME1223



**Extension set microbore
bifuse 2 female luers spin
male luer lock. Not made with
DEHP. L: 4 in L: 10 cm PV:
0.2 mL content sterile**

SKU: ME1222




**Extension set roller clamp,
spin male luer lock. Not made
with DEHP. L: 30 in L: 76 cm
PV: 5 mL. content sterile**

SKU: ME1217




**Extension set microbore spin
male luer lock. Not made with
DEHP. L: 72 in L: 183 cm PV:
0.4 mL content sterile**

EXT SET MICROBORE 84" PINCH
CLMP FIX LL
SKU: ME2059




**Extension set minibore
removable slide clamp(s) fixed
male luer lock. Not made with
DEHP. L: 60 in L: 152 cm PV:
1.6 mL content sterile**


SKU: ME2050




**Extension set minibore
removable slide clamp(s)
standard injection t-connector
fixed male luer lock. Not made
with DEHP. L: 6 in L: 15 cm
PV: 0.3 mL content sterile**
SKU: ME1214




**Extension set 1 clear needle-
free access connector
removable slide clamp(s) spin
male luer lock. Not made with
DEHP. Content sterile**
SKU: ME1213



**Extension set minibore
removable slide clamp(s) fixed
male luer lock. Not made with
DEHP. L: 36 in L: 91 cm PV:
1.3 mL content sterile**
SKU: ME1110




**Extension set slide clamp(s)
spin male luer lock. Not made
with DEHP. L: 20 in L: 51 cm
PV: 2.4 mL content sterile**
SKU: ME1069




**Extension set slide clamp(s)
spin male luer lock. Not made
with DEHP. L: 30 in L: 76 cm
PV: 4.9 mL content sterile**
SKU: ME1055




**Extension set minibore
removable slide clamp(s) spin
male luer lock. Not made with
DEHP. L: 12 in L: 30.5 cm
PV: 0.4 mL content sterile**
SKU: ME1048



**Extension set minibore
removable slide clamp(s)
standard injection t-connector
male luer slip adapter. Not
made with DEHP. L: 4 in L:
10.2 cm PV: 0.33 mL content
sterile**
SKU: ME1046




**Extension set minibore
removable slide clamp(s) spin
male luer lock. Not made with
DEHP. L: 6 in L: 15 cm PV:
0.3 mL content sterile**
SKU: ME1032



Extension set minibore trifuse
4 removable slide clamp(s) 3
female luers spin male luer
lock. Not made with DEHP. L:
6 in L: 15 cm PV: 0.5 mL
content sterile

SKU: ME1021



Extension set minibore bifuse
3 removable slide clamp(s) 2
female luers spin male luer
lock. Not made with DEHP. L:
6 in L: 15 cm PV: 0.6 mL
content sterile

SKU: ME1001



Extension set 4-way stopcock
(off indicator) spin male luer
lock. Not made with DEHP. L:
36 in L: 91 cm PV: 5.5 mL
content sterile

SKU: M4058




Extension set 3-way stopcock
(off indicator) spin male luer
lock. Not made with DEHP. L:
20 in L: 51 cm PV: 3.2 mL
content sterile

SKU: M4026



**Extension set smallbore spin
male luer lock. Not made with
DEHP. L: 63 in L: 160 cm PV:
1.6 mL fluid path sterile**

MIC EXT SET
SKU: C30006



**SmartSite™ Extension Set,
smallbore, 1 removable
SmartSite™ Needle-Free
Connector, slide clamps,
SmartSite™ T-Connector spin
male luer lock**

SKU: C25013E



**Extension set smallbore anti-
siphon valve 2 male luers
fixed male luer lock. Not made
with DEHP. L: 5 in L: 13 cm
PV: 0.5 mL fluid path sterile**

DEHP FREE 4 IN EXT SET
W/ANTSPHN VLV DET
SKU: C25012



**Extension set microbore pinch
clamp(s) fixed male luer lock.
Not made with DEHP. L: 84 in
L: 213 cm PV: 0.7 mL fluid
path sterile**

SKU: C25006



**Extension set pinch clamp(s)
2-piece male luer lock. L: 8 in
L: 20 cm PV: 1 mL fluid path
sterile**

SKU: C20025



**Extension set pinch clamp(s)
2-piece male luer lock. Not
made with DEHP. L: 20 in L:
51 cm PV: 3 mL fluid path
sterile**

SKU: C20022




**Extension set smallbore fixed
male luer lock. Not made with
DEHP. L: 30 in L: 76 cm PV:
2 mL fluid path sterile**

EXT SET SMALLBORE 40" SLIDE
CLMP LUER SLIP
SKU: C20016



**Extension set low sorbing 2-
piece male luer lock. Not
made with DEHP. L: 30 in L:
76 cm PV: 4 mL fluid path
sterile**

SKU: C20014



**BD Q-Syte™ Needle-free
Connector, 50 cm (20 in.) "2y"
macro bore extension set and
spin nut**

SKU: 385166



**SmartSite™ Extension Set, 2
pinch clamps, 2
SmartSite™ Needle-free
Valves, 7", 32" from fixed
male luer lock**

SKU: 37262E



**Extension set microbore slide
clamp(s) fixed male luer lock.
Not made with DEHP. L: 60 in
L: 152 cm PV: 0.55 mL fluid
path sterile**

SKU: 30914



**Extension set microbore slide
clamp(s) fixed male luer lock.
Not made with DEHP. L: 36 in
L: 91 cm PV: 0.4 mL fluid
path sterile**


SKU: 30910



Extension set microbore anti-siphon valve pinch clamp(s) fixed male luer lock. Not made with DEHP. L: 62 in L: 156 cm PV: 1mL fluid path sterile
EXT DEHP W/ANTISIPHON VLV
SKU: 30862-0500



Extension set microbore pinch clamp(s) fixed male luer lock. Not made with DEHP. L: 60 in L: 152 cm PV: 0.5 mL fluid path sterile
SKU: 30822



SmartSite™ Extension Set, smallbore, pinch clamps, 1 SmartSite™ Needle-Free Valve, 9" from spin male luer lock
SKU: 30212E



SmartSite™ Extension Set, pinch clamps, 1 SmartSite™ Needle-Free Valve, 9" from spin male luer lock
SKU: 30202E



**SmartSite™ Extension Set, 3
ganged 4-way stopcocks, 1
SmartSite™ Needle-Free
Valve, 6" from spin male luer
lock**

SKU: 25208E



**SmartSite™ Extension Set,
pinch clamps, 2 ganged 4-way
stopcocks, 1 SmartSite™
Needle-Free Valve, 8" from
spin male luer lock**

SKU: 25207E



**Extension set 4-way stopcock
spin male luer lock. Not made
with DEHP. L: 32 in L: 81 cm
PV: 4 mL fluid path sterile**

SKU: 25206




**Extension set 3-way stopcock
spin male luer lock. Not made
with DEHP. L: 32 in L: 81 cm
PV: 4 mL fluid path sterile**

SKU: 25205



**Extension set 3-way stopcock
spin male luer lock. Not made
with DEHP. L: 20 in L: 51 cm
PV: 3 mL fluid path sterile**

SKU: 25203




**SmartSite™ Extension Set,
pinch clamps, 4-way stopcock
(off indicator), 1
SmartSite™ Needle-Free
Valve, 6" from 2-piece male
luer lock**

SKU: 25202E



**SmartSite™ Extension Set,
pinch clamps, 1
SmartSite™ Needle-Free
Valve, 10" from spin male luer
lock**

SKU: C25004E



**SmartSite™ Extension Set,
smallbore, 4
SmartSite™ Needle-Free
Valves, 6", 14", 218", 226"
from fixed male luer lock**

SKU: C25002E



SmartSite™ Extension Set, 1 removable SmartSite™ Needle-Free Connector, slide clamps, spin male luer lock

SKU: 22059E




SmartSite™ Extension Set, microbore 1 SmartSite™ Needle-Free connector, slide clamp

SKU: 22005E-07T



SmartSite™ Extension Set, pinch clamp, 3 port manifold (3 one way valves), 1 SmartSite™ Needle-Free Valves, 6" from 2-piece male luer lock

SKU: 20519E



SmartSite™ Extension Set, pinch clamps, 4-way stopcock, 2 one way valves, 2 male luers

SKU: 20511E




**SmartSite™ Extension Set,
low sorbing, 0.2 micron filter
(1.72 sq. in low protein
binding), slide clamps**

SKU: 20350ET




**SmartSite™ Extension Set,
low sorbing, 0.2 micron filter
(1.72 sq. in low protein
binding), pinch clamps**

SKU: 20350E




**SmartSite™ Extension Set,
smallbore, trifuse 1 removable
SmartSite™ Needle-Free
Connector, 4 slide clamps, 2
female luers spin male luer
lock**

SKU: 20338E



**SmartSite™ Extension Set,
bifuse 2 SmartSite™ Needle-
Free Connectors, 2 slide
clamps, spin male luer lock**

SKU: 20159E



**SmartSite™ Extension Set,
smallbore, check valve trifuse
3 SmartSite™ Needle-Free
Connectors, 3 slide clamps,
spin male luer lock**


Extension set smallbore check valve
trifuse 3 SmartSite™ needle-free
connector(s) 3 slide clamp(s) spin male
luer lock. L: 7.5 in L: 19 cm PV: 1.7 mL
fluid path sterile.

SKU: 20062E



**SmartSite™ Extension Set, 1
SmartSite™ Needle-Free
Connector, slide clamps, spin
male luer lock**

SKU: 20059E




**SmartSite™ Extension Set,
smallbore slide clamps,
SmartSite™ T-Connector
male luer slip adapter**

SKU: 20051E




BD SmartSite™ Bag 500 mL

SKU: 200500




**SmartSite™ Extension Set,
smallbore, 1
SmartSite™ Needle-Free
Connector, slide clamps, spin
male luer lock**

SKU: 20049E




**SmartSite™ Extension Set,
smallbore slide clamps,
SmartSite™ T-connector spin
male luer lock**

SKU: 20043E




**SmartSite™ Extension Set,
pressure rated smallbore slide
clamps SmartSite™ T-
connector spin male luer lock**

SKU: 20041E



**SmartSite™ Extension Set,
pressure rated smallbore 1
SmartSite™ Needle-Free
Connector, slide clamps, spin
male luer lock**

SKU: 20039E



**SmartSite™ Extension Set,
smallbore trifuse 3
SmartSite™ Needle-Free
Connectors, 3 slide clamps,
fixed male luer lock**

SKU: 20038E



**SmartSite™ Extension Set,
smallbore 0.2 micron small
volume filter (0.25 sq. in low
protein binding), slide clamps**

SKU: 20029E




**SmartSite™ Extension Set, 0.2
micron filter (1.72 sq. in low
protein binding), pinch
clamps**

SKU: 20027E



BD SmartSite™ Bag 250 mL

SKU: 200250



**SmartSite™ Extension Set
smallbore, 1
SmartSite™ Needle-Free
Connector, slide clamps, spin
male luer lock J-loop**

SKU: 20021E



**SmartSite™ Extension Set
pressure rated smallbore
bifuse, 2 SmartSite™ Needle-
Free Connectors, 2 slide
clamps, spin male luer lock**

SKU: 20019E



BD SmartSite™ Bag 100 mL

SKU: 200100




**Extension set 2-piece male
luer lock. Not made with
DEHP. L: 21 in L: 53 cm PV:
2.9 mL fluid path sterile**

SKU: 11608375



**SmartSite™ Extension Set,
pinch clamps, 1
SmartSite™ Needle-Free
Valves, 4" from spin male luer
lock**

SKU: 11590099



**Extension set smallbore slide
clamp(s) spin male luer lock.
Not made with DEHP. L: 62 in
L: 158 cm PV: 1.5 mL fluid
path sterile**

EXT SET DEHP FREE 158


SKU: 11088483



**Extension set smallbore pinch
clamp(s) spin male luer lock.
Not made with DEHP. L: 110
in L: 279 cm PV: 6 mL fluid
path sterile**

EXT SET LIGHT RES TBG

SKU: 10896354




**Extension set smallbore
trifuse 4 slide clamp(s) 3
female luers spin male luer
lock. Not made with DEHP. L:
7 in L: 18 cm PV: 0.5 mL fluid
path sterile**

SKU: 10839681



**Extension set smallbore bifuse
2 slide clamp(s) 2 female luers
spin male luer lock. Not made
with DEHP. L: 7 in L: 18 cm
PV: 0.3 mL fluid path sterile**
SKU: 10794983



**Extension set smallbore slide
clamp(s) spin male luer lock.
Not made with DEHP. L: 6 in
L: 15 cm PV: 0.3 mL fluid
path sterile**
SKU: 10373772



**Extension set microbore
amber slide clamp(s) fixed
male luer lock. Not made with
DEHP. L: 62 in L: 156 cm PV:
0.6 mL fluid path sterile**
EXT SET MICRO DEHP-FREE
SKU: 10241342




**Extension set smallbore pinch
clamp(s) fixed male luer lock.
Not made with DEHP. L: 234
in L: 594 cm PV: 12 mL fluid
path sterile**
EXT SET MICRO DEHP-FREE
SKU: 10031776



**Extension set slide clamp(s) 2-
piece male luer lock. L: 8 in
L: 20 cm PV: 0.9 mL fluid
path sterile**

SKU: 10026685



**Extension set smallbore
quadfuse 4 slide clamp(s) 4
female luers spin male luer
lock. Not made with DEHP. L:
8 in L: 20 cm PV: 0.7 mL fluid
path sterile**

SKU: 10015817



**Extension set smallbore 1.2
micron small volume filter
(0.25 sq. in low protein
binding) slide clamp(s) fixed
male luer lock. Not made with
DEHP. L: 7 in L: 18 cm PV:
0.5 mL fluid path sterile**

SKU: 20129E



**SmartSite™ Extension Set, 1.2
micron filter (1.72 sq. in low
protein binding), pinch
clamps, 1 SmartSite™ Needle-
Free Valve**

SKU: 20127E




**SmartSite™ Extension Set
smallbore bifuse, 2 removable
SmartSite™ Needle-Free
Connectors, 2 slide clamps,
spin male luer lock**

SKU: 10015488




**Extension set pinch clamp(s)
2-piece male luer lock. Not
made with DEHP. L: 8 in L:
20 cm PV: 1 mL fluid path
sterile**

SKU: 10015485



**Extension set smallbore slide
clamp(s) VersaSafe™ t-
connector male luer slip
adapter. Not made with
DEHP. L: 6 in L: 15 cm PV:
0.2 mL fluid path sterile**

SKU: 10014813



**SmartSite™ Extension Set,
trifuse, 2 SmartSite™ Needle-
Free Connectors, 4 slide
clamps, 2-piece fixed male
luer lock**

SKU: 10014692



Extension set slide clamp(s) 4-way stopcock 2 one way valves 2-piece male luer lock. Not made with DEHP. L: 45 in L: 114 cm PV: 7.5 mL fluid path sterile

SKU: 10014472



SmartSite™ Extension Set, low sorbing 1 SmartSite™ Needle-Free Connector 0.2 micron filter (1.72 sq. in low protein binding), pinch clamp

SKU: 10013902




Extension set 1.2 micron filter (1.72 sq. in low protein binding) spin male luer lock. Not made with DEHP. L: 17 in L: 43 cm PV: 4 mL fluid path sterile

SKU: 10012866




Extension set pinch clamp(s) 2-piece male luer lock. L: 60 in L: 152 cm PV: 7 mL fluid path sterile

SKU: 10012213




**SmartSite™ Extension Set
bifuse, 1 removable
SmartSite™ Needle-Free
Connector, 2 slide clamps,
spin male luer lock**

SKU: 10012143



**Extension set smallbore 0.2
micron small volume filter
(0.25 sq. in low protein
binding) slide clamp(s) spin
male luer lock. Not made with
DEHP. L: 7 in L: 18 cm PV:
0.5 mL fluid path sterile**

SKU: 10011865



**SmartSite™ Extension Set
pressure rated smallbore, 1
removable SmartSite™
Needle-Free Connector, slide
clamp, spin male luer lock**

SKU: 10011253




**SmartSite™ Extension Set,
smallbore, 1 removable
SmartSite™ Needle-Free
Connector, slide clamps, fixed
male luer lock**


SKU: 10010983




**SmartSite™ Extension Set, 1
SmartSite™ Needle-free
Connector slide clamp, spin
male luer lock, L: 9 in L: 23
cm, Priming Volume: 1.3 mL**
SKU: 10010912



**Extension set microbore 0.2
micron small volume filter
(0.25 sq. in low protein
binding) slide clamp(s) fixed
male luer lock. DEHP is not
part of the material
formulation. L: 60 in L: 152
cm PV: 1 mL fluid path sterile**
SKU: 10010570



**SmartSite™ Extension Set,
smallbore, 1
SmartSite™ Needle-Free
Connector, slide clamp, fixed
male luer lock**
SKU: 10010511



Gravity administration set

SKU: MZ8002



Bag access spike with cap, clear needle-free access connector

SKU: MY3003



Bag access spike with cap, back check valve, clear needle-free access connector

SKU: MY3001




Gravity administration set

SKU: MFS124




BD Alaris™ LVP 20D DEHP AMBER SS 0.2M

AS LVP 20D DEHP AMBER SS 0.2M
SKU: C24101E



Monoject™ SYRINGE BARREL CAP

SKU: 8881682101



Cardinal Health™ Monoject™ SYRINGE BARREL ASSEMBLY

SKU: 8881135609



Gravity administration set

SKU: 42433E



Gravity administration set

SKU: 42293E



Gravity administration set

SKU: 42233E



Gravity set

SKU: 42000



Gravity administration set

SKU: 41177E



Gravity administration set

SKU: 41134E



BD Alaris™ PCA ASV MICROBORE

SKU: 30883



BD Alaris™ PCA ASV MICROBORE CV

PCA ASV MICROBORE CV

SKU: 30873

BD Alaris™ PCA SS CV


PCA SS CV

SKU: 30843E




Alaris™ MS III HALF DEHP SS

SKU: 28117E



Alaris™ MS III BLD200M 20D DEHP SS

SKU: 28080E



Alaris MedSystem III™ infusion pump administration set

SKU: 28034



BD Alaris™ LVP BLD180M 15D SS

AS LVP BLD180M 15D DEHP 1SS

SKU: 2477-0007



BD Alaris™ Pump Module administration set

SKU: 2477-0000



BD Alaris™ LVP SS BAG 20D LOW SORB

AS LVP SS BAG 20D LOW SORB

SKU: 2466-0007



BD Alaris™ LVP SS BAG 20D LOW SORB SS 0.2M

AS LVP SS BAG 20D LOW SORB SS

0.2M

SKU: 2465-0007



BD Alaris™ LVP SS BAG 20D LOW SORB 2SS TEX CV

AS LVP SS BAG 20D LOW SORB 2SS

TEX CV

SKU: 24601-B007T



BD Alaris™ LVP 20D LOW SORB 3SS CV

AS LVP 20D LOW SORB 3SS CV

SKU: 24600-0007

BD Alaris™ LVP 20D DEHP SS

SKU: 2455-0500



BD Alaris™ LVP 10D 3SS 2CV

SKU: 2452-0007



BD Alaris™ LVP BUR 60D 2SS BV

AS LVP BUR 60D 2SS BV

SKU: 2447-0007



BD Alaris™ LVP BUR 60D SMBORE 3SS

AS LVP BUR 60D DEHP SMBORE 3INJ

SKU: 2441-0007

BD Alaris™ LVP 20D SS 0.2M CV

AS LVP 20D SS 0.2M CV

SKU: 2434-0007



BD Alaris™ Pump Module administration set

SKU: 2434-0004



BD Alaris™ LVP 20D 3SS 0.2M CV

AS LVP 20D DEHP SS 0.2M

SKU: 2432-0007



BD Alaris™ Pump Module administration set

SKU: 2432-0004



BD Alaris™ Pump Module administration set

SKU: 24302-0004



Alaris™ MS III 20D DEHP 3SS 2CVC83:E83

SKU: 28493E

Alaris™ MS III 60D DEHP 2INJ 0.2M 2CV

SKU: 28273E



BD Alaris™ LVP SS BAG 20D LOW SORB SS TEX 0.2M

AS LVP SS BAG 20D LOW SORB SS
TEX 0.2M
SKU: 24301-0007T



BD Alaris™ LVP 20D DEHP 3SS CV

SKU: 2426-0500




BD Alaris™ LVP 20D 3SS CV

AS LVP 20D 3SS CV
SKU: 2426-0007



BD Alaris™ Pump Module administration set

SKU: 2426-0004



BD Alaris™ LVP 10D 4SS 4SS 2G-4WSPK CV

AS LVP 10D DEHP 4SS 4WSPK CV

SKU: 2423-0007



BD Alaris™ LVP 20D 2 SS

SKU: 24201-0007



BD Alaris™ LVP 20D DEHP 2SS CV

SKU: 2420-0500



BD Alaris™ LVP 20D 2SS CV

AS LVP 20D 2VS CV

SKU: 2420-0007



BD Alaris™ Pump Module administration set

SKU: 2420-0004

BD Alaris™ LVP HALF SS

AS LVP HALF SS

SKU: 2403-0007



BD Alaris™ LVP 20D 3SS TEX CV

AS LVP 20D 3SS TEX CV

SKU: 24010-0007T



BD Alaris™ Pump Module administration set

SKU: 2401-0004



BD Alaris™ Pump Module administration set


SKU: 24009-0007T



BD Alaris™ LVP 20D SS CV

AS LVP 20D SS CV

SKU: 24001-0007



BD Alaris™ Pump Module administration set

SKU: 2280-0000

BD Alaris™ LVP BLD200M 20D DEHP

SKU: 2278-0500



BD Alaris™ LVP SS BAG 20D LOW SORB TEX

AS LVP SS BAG 20D TEX
SKU: 22602-B007T



BD Alaris™ LVP 20D LOW SORB TEX

AS LVP 20D LOW SORB TEX 47
SKU: 22600-0007T



BD Alaris™ LVP 20D LOW SORB

AS LVP 20D LOW SORB
SKU: 2260-0500



BD Alaris™ LVP 20D YELLOW MICROBORE

SKU: 2206-0007



BD Alaris™ LVP 20D

AS LVP 20D

SKU: 2204-0007



BD Alaris™ LVP 20D DEHP

AMBER 15M


SKU: 2203-0500



BD Alaris™ LVP 20D 1.2M

AS LVP 20D 1.2M

SKU: 2202-0007



BD Alaris™ LVP SS BAG 20D

TEX

Alaris Medley Infusion Set PVC/nDEHP

macro bore tubing General spike

Integrated spike w DC SS Bag access

20DP 1RC Texium 115cm

SKU: 22000-B007T



BD Alaris™ Pump Module administration set

SKU: 2200-0500

BD Alaris™ Pump Module administration set

SKU: 2126-0007



BD Alaris™ LVP 20D DEHP
2VS CV
SKU: 2120-0500



BD Alaris™ Pump Module
administration set
SKU: 2120-0007

100

NDC 17271-750-07

1,000 mL

freeflex®

Sterile Water for Injection, USP

DO NOT GIVE INTRAVENOUSLY.
FOR DRUG DILUENT USE ONLY.

100

For intravenous use.
Contains no antimicrobial or other substance. Additives should be used promptly with adequate precautions for maintaining sterility, pH 5.5 (5.0 to 7.0).

300

WARNINGS: NOT ISOTONIC. HEMOLYTIC.
Do not give intravenously unless rendered nearly isotonic.
Single-Dose Container. Discard Unused Portion.

400

Additives may be incompatible. Consultation with pharmacist is recommended. When introducing additives, use aseptic technique, mix thoroughly and do not store.
Usual dosage: See package insert.
The overwrap is a moisture barrier.
Use immediately once removed from overwrap.

500

Store at 20° to 25°C (68° to 77°F) [see USP Controlled Room Temperature]. Protect from freezing. Sterile, nonpyrogenic. Use only if solution is clear and container is undamaged. Must not be used in series connections.

600

The container closure is not made with natural rubber latex.
Non-PVC, Non-DEHP, Sterile.

700


(01)00317271750078

800

 **BD**
Becton, Dickinson and Company
1 Becton Drive
Franklin Lakes, NJ 07417 USA
EXP
For product inquiry: 1-800-523-0502
Distributed by BD. Manufactured by Fresenius Kabi.
Made in Germany.
403572

900

1727175007
074409100 US BD

Sterile Water Injection, USP
1000ml
NDC 17271-750-07
SKU: 1727175007

100

NDC 17271-721-07

1,000 mL

freeflex®

10% Dextrose Injection, USP

100 grams per 1,000 mL
(100 mg per mL)

200

For intravenous use.
Each 100 mL contains: Dextrose monohydrate 10 g; water for injection, 100 mL.

300

Hypertonic.
505 mOsmol/L (CALC.) pH 4.5 (3.2 to 6.5).
Single-dose Container. Discard Unused Portion.
Additives may be incompatible. Consult with pharmacist. When introducing additives, use aseptic technique, mix thoroughly and do not store. Do not administer simultaneously with blood.
Usual dosage: See package insert.
The overwrap is a moisture barrier.
Use immediately once removed from overwrap.

400

STORE AT: 20° to 25°C (68° to 77°F) [see USP Controlled Room Temperature]. Protect from freezing. Avoid excessive heat.
The container closure is not made with natural rubber latex.
Non-PVC, Non-DEHP, Sterile.

500


(01)00317271721078

600

 **BD**
Becton, Dickinson and Company
1 Becton Drive
Franklin Lakes, NJ 07417 USA
EXP
For product inquiry: 1-800-523-0502
Distributed by BD. Manufactured by Fresenius Kabi.
Made in Germany.
403566

700

1727172107

10% Dextrose Injection, USP
1000ml
NDC 17271-721-07
SKU: 1727172107

100

NDC 17271-721-06

500 mL

freeflex®

10% Dextrose Injection, USP

50 grams per 500 mL
(100 mg per mL)

200

For intravenous use.
Each 100 mL contains: Dextrose monohydrate 10 g; water for injection, 100 mL.

300

Hypertonic.
505 mOsmol/L (CALC.) pH 4.5 (3.2 to 6.5).
Single-dose Container. Discard Unused Portion.
Additives may be incompatible. Consult with pharmacist. When introducing additives, use aseptic technique, mix thoroughly and do not store. Do not administer simultaneously with blood.
Usual dosage: See package insert.
The overwrap is a moisture barrier.
Use immediately once removed from overwrap.

400

STORE AT: 20° to 25°C (68° to 77°F) [see USP Controlled Room Temperature]. Protect from freezing. Avoid excessive heat.
The container closure is not made with natural rubber latex.
Non-PVC, Non-DEHP, Sterile.

500


(01)00317271721061

600

 **BD**
Becton, Dickinson and Company
1 Becton Drive
Franklin Lakes, NJ 07417 USA
EXP
For product inquiry: 1-800-523-0502
Distributed by BD. Manufactured by Fresenius Kabi.
Made in Germany.
403565

700

1727172106

Sterile Water Injection, USP
1000ml
NDC 17271-750-07
SKU: 1727175007

100

NDC 17271-721-05

250 mL

freeflex®

10% Dextrose Injection, USP

25 grams per 250 mL
(100 mg per mL)

200

For intravenous use.
Each 100 mL contains: Dextrose monohydrate 10 g; water for injection, 100 mL.

300

Hypertonic.
505 mOsmol/L (CALC.) pH 4.5 (3.2 to 6.5).
Single-dose Container. Discard Unused Portion.
Additives may be incompatible. Consult with pharmacist. When introducing additives, use aseptic technique, mix thoroughly and do not store. Do not administer simultaneously with blood.
Usual dosage: See package insert.
The overwrap is a moisture barrier.
Use immediately once removed from overwrap.

400

STORE AT: 20° to 25°C (68° to 77°F) [see USP Controlled Room Temperature]. Protect from freezing. Avoid excessive heat.
The container closure is not made with natural rubber latex.
Non-PVC, Non-DEHP, Sterile.

500


(01)00317271721054

600

 **BD**
Becton, Dickinson and Company
1 Becton Drive
Franklin Lakes, NJ 07417 USA
EXP
For product inquiry: 1-800-523-0502
Distributed by BD. Manufactured by Fresenius Kabi.
Made in Germany.
403564

700

1727172105

10% Dextrose Injection, USP
1000ml
NDC 17271-721-07
SKU: 1727172107

10% Dextrose Injection, USP
500ml

NDC 17271-721-06
SKU: 1727172106

100

5% Dextrose Injection, USP

200

50 grams per 1,000 mL
(50 mg per mL)

300

400

500

600

700

800

900

For intravenous use. Rx only

Each 100 mL contains: Dextrose monohydrate, 5 g; water for injection, 100 mL.

252 mOsmol/LITER (CALC.) pH 4.3 (3.2 to 6.5).

Single Dose Only. Discard Unused Portion.

Dextrose solutions without salts should not be used in blood transfusions because of possible Rouleau formation. Additives may be incompatible. Consult with pharmacist. When introducing additives, use aseptic technique, mix thoroughly and do not store.

Usual dosage: See package insert. The overwrap is a moisture barrier. Use immediately once removed from overwrap.

STORE AT: 20° to 25°C (68° to 77°F) [see USP Controlled Room Temperature]. Protect from freezing.

The container closure is not made with natural rubber latex. Non-PVC, Non-DEHP, Sterile.



(01) 0 03 1727172007 1 NDC: 17271-720-07

BD

LOT

Becton, Dickinson and Company
1 Becton Drive
Franklin Lakes, NJ 07417 USA
For product inquiry: 1-800-523-0502
Distributed by BD. Manufactured by Fresenius Kabi.
Made in Germany.

EXP

403239
1727172007
0743341/00 US

10% Dextrose Injection, USP
250ml

NDC 17271-721-05
SKU: 1727172105

100

5% Dextrose Injection, USP

200

25 grams per 500 mL
(50 mg per mL)

300

400

500

600

700

800

900

For intravenous use. Rx only

Each 100 mL contains: Dextrose monohydrate, 5 g; water for injection, 100 mL.

252 mOsmol/LITER (CALC.) pH 4.3 (3.2 to 6.5).

Single Dose Only. Discard Unused Portion.

Dextrose solutions without salts should not be used in blood transfusions because of possible Rouleau formation. Additives may be incompatible. Consult with pharmacist. When introducing additives, use aseptic technique, mix thoroughly and do not store.

Usual dosage: See package insert. The overwrap is a moisture barrier. Use immediately once removed from overwrap.

STORE AT: 20° to 25°C (68° to 77°F) [see USP Controlled Room Temperature]. Protect from freezing.

The container closure is not made with natural rubber latex. Non-PVC, Non-DEHP, Sterile.



(01) 0 03 1727172006 4 NDC: 17271-720-06

BD

LOT

Becton, Dickinson and Company
1 Becton Drive
Franklin Lakes, NJ 07417 USA
For product inquiry: 1-800-523-0502
Distributed by BD. Manufactured by Fresenius Kabi.
Made in Germany.

EXP

403237
1727172006
0743361/00 US

5% Dextrose Injection, USP
1000mL

NDC 17271-720-07
SKU: 1727172007



5% Dextrose Injection, USP
500mL

NDC 17271-720-06
SKU: 1727172006



BD Alaris™ Pump Module
administration set

SKU: 22604-B007

BD Alaris™ LVP SS BAG 20D
LOW SORB TEX

AS LVP SS BAG 20D LOW SORB TEX
47

SKU: 22603-B007T

5% Dextrose Injection, USP

12.5 grams per 250 mL
(50 mg per mL)

For intravenous use.
Rx only

Each 100 mL contains: Dextrose monohydrate, 5 g; water for injection, 100 mL. 252 mOsmol/LITER (CALC.) pH 4.3 (3.2 to 6.5). Single Dose Only. Discard Unused Portion. Dextrose solutions without salts should not be used in blood transfusions because of possible Rouleau formation. Additives may be incompatible. Consult with pharmacist. When introducing additives, use aseptic technique, mix thoroughly and do not store. Usual dosage: See package insert. The overwrap is a moisture barrier. Use immediately once removed from overwrap. STORE AT: 20° to 25°C (68° to 77°F) [see USP Controlled Room Temperature]. Protect from freezing. The container closure is not made with natural rubber latex. Non-PVC, Non-DEHP, Sterile.

403230 1727172005
FDH 0072 01-62-04-006

BD

LOT

Becton, Dickinson and Company
1 Becton Drive
Franklin Lakes, NJ 07417 USA
For product inquiry: 1-800-523-0502
Distributed by BD. Manufactured by Fresenius Kabi. Made in Norway.

EXP

5% Dextrose Injection, USP
250mL
NDC 17271-720-05
SKU: 1727172005

5% Dextrose Injection, USP

2.5 grams per 50 mL
(50 mg per mL)

For intravenous use.
Rx only

Each 100 mL contains: Dextrose monohydrate, 5 g; water for injection, 100 mL. 252 mOsmol/LITER (CALC.) pH 4.3 (3.2 to 6.5). Single Dose Only. Discard Unused Portion. Dextrose solutions without salts should not be used in blood transfusions because of possible Rouleau formation. Additive compatibility, consult pharmacist. When introducing additives, use aseptic technique, mix thoroughly and do not store. Usual dosage: See package insert. The overwrap is a moisture barrier. Use immediately once removed from overwrap. STORE AT: 20° to 25°C (68° to 77°F) [see USP Controlled Room Temperature]. Protect from freezing. The container closure is not made with natural rubber latex. Non-PVC, Non-DEHP, Sterile.

403562 FDE 0076 01-69-04-006 1 2 3 4 5 6 7 8 9 0

BD

LOT

Becton, Dickinson and Company
1 Becton Drive
Franklin Lakes, NJ 07417 USA
For product inquiry: 1-800-523-0502
Distributed by BD. Manufactured by Fresenius Kabi. Made in Norway.

EXP

5% Dextrose Injection, USP
50ml
NDC 17271-720-02
SKU: 1727172002

freeflex[®] NDC 17271-710-06 500 mL

Lactated Ringer's Injection, USP

For intravenous use. Rx only

Each 100 mL contains:
Sodium Chloride, USP 600 mg Na⁺ 130
Sodium Lactate, USP 310 mg K⁺ 4
Potassium Chloride, USP 30 mg Ca²⁺ 2.7
Calcium Chloride, USP 20 mg Cl⁻ 109
273 mOsmol/L (calc.) pH 6.5 (6.0 to 7.5)
Single Dose Only. Discard Unused Portion.

Not for use in the treatment of lactic acidosis.

Additives may be incompatible. Consult with pharmacist. When introducing additives, use aseptic technique, mix thoroughly and do not store. Use only if solution is clear and container is undamaged. Must not be used in series connections. Do not administer simultaneously with blood.

Usual dosage: See package insert. The overwrap is a moisture barrier. Use immediately once removed from overwrap.

STORE AT: 20° to 25°C (68° to 77°F) [see USP Controlled Room Temperature]. Avoid excessive heat.

The container closure is not made with natural rubber latex. Non-PVC, Non-DEHP, Sterile.

1234567890

BD

LOT

Becton, Dickinson and Company
1 Becton Drive
Franklin Lakes, NJ 07417 USA
For product inquiry: 1-800-523-0502
Distributed by BD. Manufactured by Fresenius Kabi. Made in Norway.

EXP

5% Dextrose Injection, USP

5 grams per 100 mL
(50 mg per mL)

For intravenous use.
Rx only

Each 100 mL contains: Dextrose monohydrate, 5 g; water for injection, 100 mL. 252 mOsmol/LITER (CALC.) pH 4.3 (3.2 to 6.5). Single Dose Only. Discard Unused Portion. Dextrose solutions without salts should not be used in blood transfusions because of possible Rouleau formation. Additive compatibility, consult pharmacist. When introducing additives, use aseptic technique, mix thoroughly and do not store. Usual dosage: See package insert. The overwrap is a moisture barrier. Use immediately once removed from overwrap. STORE AT: 20° to 25°C (68° to 77°F) [see USP Controlled Room Temperature]. Protect from freezing. The container closure is not made with natural rubber latex. Non-PVC, Non-DEHP, Sterile.

403225 1727172003 FDH 0075 01-62-04-005

BD

LOT

Becton, Dickinson and Company
1 Becton Drive
Franklin Lakes, NJ 07417 USA
For product inquiry: 1-800-523-0502
Distributed by BD. Manufactured by Fresenius Kabi. Made in Norway.

EXP

5% Dextrose Injection, USP
100mL
NDC 17271-720-03
SKU: 1727172003

freeflex[®] NDC 17271-710-07 1,000 mL

Lactated Ringer's Injection, USP

For intravenous use. Rx only

Each 100 mL contains:
Sodium Chloride, USP 600 mg Na⁺ 130
Sodium Lactate, USP 310 mg K⁺ 4
Potassium Chloride, USP 30 mg Ca²⁺ 2.7
Calcium Chloride, USP 20 mg Cl⁻ 109
273 mOsmol/L (calc.) pH 6.5 (6.0 to 7.5)
Single Dose Only. Discard Unused Portion.

Not for use in the treatment of lactic acidosis.

Additives may be incompatible. Consult with pharmacist. When introducing additives, use aseptic technique, mix thoroughly and do not store. Use only if solution is clear and container is undamaged. Must not be used in series connections. Do not administer simultaneously with blood.

Usual dosage: See package insert. The overwrap is a moisture barrier. Use immediately once removed from overwrap.

STORE AT: 20° to 25°C (68° to 77°F) [see USP Controlled Room Temperature]. Avoid excessive heat.

The container closure is not made with natural rubber latex. Non-PVC, Non-DEHP, Sterile.

403509 FDH 2344 01-62-12-032

BD

LOT

Becton, Dickinson and Company
1 Becton Drive
Franklin Lakes, NJ 07417 USA
For product inquiry: 1-800-523-0502
Distributed by BD. Manufactured by Fresenius Kabi. Made in Norway.

EXP

Lactated Ringer's Injection, USP 1000mL
NDC 17271-710-07
SKU: 1727171007

freeflex[®] NDC 17271-710-05 250 mL

Lactated Ringer's Injection, USP

For intravenous use. Rx only

Each 100 mL contains:
Sodium Chloride, USP 600 mg Na⁺ 130
Sodium Lactate, USP 310 mg K⁺ 4
Potassium Chloride, USP 30 mg Ca²⁺ 2.7
Calcium Chloride, USP 20 mg Cl⁻ 109
273 mOsmol/L (calc.) pH 6.5 (6.0 to 7.5)
Single Dose Only. Discard Unused Portion.

Not for use in the treatment of lactic acidosis.

Additives may be incompatible. Consult with pharmacist. When introducing additives, use aseptic technique, mix thoroughly and do not store. Use only if solution is clear and container is undamaged. Must not be used in series connections. Do not administer simultaneously with blood.

Usual dosage: See package insert. The overwrap is a moisture barrier. Use immediately once removed from overwrap.

STORE AT: 20° to 25°C (68° to 77°F) [see USP Controlled Room Temperature]. Avoid excessive heat.

The container closure is not made with natural rubber latex. Non-PVC, Non-DEHP, Sterile.

1234567890

BD

LOT

Becton, Dickinson and Company
1 Becton Drive
Franklin Lakes, NJ 07417 USA
For product inquiry: 1-800-523-0502
Distributed by BD. Manufactured by Fresenius Kabi. Made in Norway.

EXP

Lactated Ringer's Injection,
USP 500mL

NDC 17271-710-06
SKU: 1727171006

100

200

300

400

freeflex

500 mL

3% Sodium Chloride Injection, USP

15 grams per 500 mL
(30 mg per mL)

Hypertonic Saline

For intravenous use.

Each 100 mL contains: Sodium Chloride, USP 3 g in water for injection. pH may have been adjusted with hydrochloric acid.

Sodium Chloride 513 mEq/L

1027 mOsmol/L (CALC.) pH 5.0 (4.5 to 7.0)

For intravenous use.

Single-Dose Container. Discard Unused Portion.

Additives may be incompatible. Consult with pharmacist, if available.

When introducing additives, use aseptic technique, mix thoroughly and do not store.

Usual dosage: See package insert.

Store unit in moisture barrier overwrap.

Use immediately once removed from overwrap.

Store at 20° to 25°C (68° to 77°F) [see USP Controlled Room Temperature] until ready to use. Avoid excessive heat.

The container closure is not made with natural rubber latex.

Non-PVC, Non-DEHP, Sterile.

BD

LOT

Becton, Dickinson and Company

EXP

1 Becton Drive

Franklin Lakes, NJ 07417 USA

For product inquiry: 1-800-523-0502

Distributed by BD. Manufactured by Fresenius Kabi.

Made in Germany

403255A

074411100 US BD

010031727170306

3% Sodium Chloride
Injection, USP 500ml

NDC 17271-703-06
SKU: 1727170306

100

200

300

400

freeflex

500 mL

0.9% Sodium Chloride Injection, USP

For intravenous use.

Rx only

Each 100 mL contains: Sodium Chloride 900 mg in water for injection.

Electrolytes per 1,000 mL:

Sodium 154 mEq

Chloride 154 mEq

308 mOsmol/LITER (CALC.) pH 4.5 to 7.0

Single Dose Only. Discard Unused Portion.

Additives may be incompatible. Consult with pharmacist. When introducing additives, use aseptic technique, mix thoroughly and do not store. Use only if solution is clear and container is undamaged. Must not be used in series connections.

Usual dosage: See package insert.

STORE AT: 20° to 25°C (68° to 77°F) [see USP Controlled Room Temperature]. Protect from freezing.

The container closure is not made with natural rubber latex.

Non-PVC, Non-DEHP, Sterile.

BD

LOT

Becton, Dickinson and Company

EXP

1 Becton Drive

Franklin Lakes, NJ 07417 USA

For product inquiry: 1-800-523-0502

Distributed by BD. Manufactured by Fresenius Kabi.

Made in Germany

403255

1727170106

0743331/00 US

010031727170106

0.9% Sodium Chloride
Injection, USP 500mL

SKU: 1727170106

Lactated Ringer's Injection,
USP 250mL

NDC 17271-710-05
SKU: 1727171005

100

200

300

400

500

600

700

freeflex

1,000 mL

0.9% Sodium Chloride Injection, USP

For intravenous use.

Rx only

Each 100 mL contains: Sodium Chloride 900 mg in water for injection.

Electrolytes per 1,000 mL:

Sodium 154 mEq

Chloride 154 mEq

308 mOsmol/LITER (CALC.) pH 4.5 to 7.0

Single Dose Only. Discard Unused Portion.

Additives may be incompatible. Consult with pharmacist. When introducing additives, use aseptic technique, mix thoroughly and do not store. Use only if solution is clear and container is undamaged. Must not be used in series connections.

Usual dosage: See package insert.

STORE AT: 20° to 25°C (68° to 77°F) [see USP Controlled Room Temperature]. Protect from freezing.

The container closure is not made with natural rubber latex.

Non-PVC, Non-DEHP, Sterile.

BD

LOT

Becton, Dickinson and Company

EXP

1 Becton Drive

Franklin Lakes, NJ 07417 USA

For product inquiry: 1-800-523-0502

Distributed by BD. Manufactured by Fresenius Kabi.

Made in Germany

403256

1727170107

0743331/00 US

010031727170107

0.9% Sodium Chloride
Injection, USP 1000mL

SKU: 1727170107

50

100

150

200

freeflex

250 mL

0.9% Sodium Chloride Injection, USP

For intravenous use.

Rx only

Each 100 mL contains: Sodium Chloride 900 mg in water for injection.

Electrolytes per 1,000 mL:

Sodium 154 mEq

Chloride 154 mEq

308 mOsmol/LITER (CALC.) pH 4.5 to 7.0

Single Dose Only. Discard Unused Portion.

Additives may be incompatible. Consult with pharmacist. When introducing additives, use aseptic technique, mix thoroughly and do not store. Use only if solution is clear and container is undamaged. Must not be used in series connections.

Usual dosage: See package insert.

STORE AT: 20° to 25°C (68° to 77°F) [see USP Controlled Room Temperature]. Protect from freezing.

The container closure is not made with natural rubber latex. Non-PVC, Non-DEHP, Sterile.

BD

LOT

Becton, Dickinson and Company

EXP

1 Becton Drive

Franklin Lakes, NJ 07417 USA

For product inquiry: 1-800-523-0502

Distributed by BD. Manufactured by Fresenius Kabi.

Made in Germany

403256

1727170107

0743331/00 US

010031727170105

0.9% Sodium Chloride
Injection, USP 250mL

SKU: 1727170105

free flex[®]

100 mL

0.9% Sodium Chloride Injection, USP

For intravenous use. Rx only

Each 100 mL contains: Sodium Chloride 900 mg in water for injection.

Electrolytes per 1,000 mL:

Sodium Chloride 154 mEq 154 mEq

308 mOsmol/LITER (CALC.) pH 4.5 to 7.0

Single Dose Only. Discard Unused Portion.

Additives may be incompatible. Consult with pharmacist. When introducing additives, use aseptic technique, mix thoroughly and do not store. Use only if solution is clear and container is undamaged. Must not be used in series connections.

Usual dosage: See package insert.

STORE AT: 20° to 25°C (68° to 77°F) [see USP Controlled Room Temperature]. Protect from freezing.

The container closure is not made with natural rubber latex. Non-PVC, Non-DEHP, Sterile.

BD

LOT

EXP

Becton, Dickinson and Company
1 Becton Drive
Franklin Lakes, NJ 07417 USA
For product inquiry: 1-800-523-0502
Distributed by BD.
Manufactured by Fresenius Kabi.
Made in Norway.

403238

1727170103

FDH 1335 01-62-19-015

(01) 0 03 1727170103 2

NDC: 17271-701-03

free flex[®]

50 mL

0.9% Sodium Chloride Injection, USP

For intravenous use. Rx only

Each 100 mL contains: Sodium Chloride 900 mg in water for injection.

Electrolytes per 1,000 mL:

Sodium Chloride 154 mEq 154 mEq

308 mOsmol/LITER (CALC.) pH 4.5 to 7.0

Single Dose Only. Discard Unused Portion.

Additives may be incompatible. Consult with pharmacist. When introducing additives, use aseptic technique, mix thoroughly and do not store. Use only if solution is clear and container is undamaged. Must not be used in series connections.

Usual dosage: See package insert.

STORE AT: 20° to 25°C (68° to 77°F) [see USP Controlled Room Temperature]. Protect from freezing.

The container closure is not made with natural rubber latex. Non-PVC, Non-DEHP, Sterile.

BD

LOT

EXP

Becton, Dickinson and Company
1 Becton Drive
Franklin Lakes, NJ 07417 USA
For product inquiry: 1-800-523-0502
Distributed by BD.
Manufactured by Fresenius Kabi.
Made in Norway.

403234

1727170102

FDH 1336 01-62-19-014

(01) 0 03 1727170102 5

NDC: 17271-701-02

0.9% Sodium Chloride Injection, USP 100 mL

SKU: 1727170103

0.9% Sodium Chloride Injection, USP 50 mL

SKU: 1727170102

free flex[®]

1,000 mL

0.45% Sodium Chloride Injection, USP

For intravenous use. Rx only

4.5 grams per 1,000 mL (4.5 mg per mL)

Each 100 mL contains: Sodium Chloride 450 mg in water for injection.

Electrolytes per 1,000 mL:

Sodium Chloride 77 mEq 77 mEq

154 mOsmol/LITER (CALC.) pH 5.6 (4.5 to 7.0)

Single-Dose Container. Discard Unused Portion.

Additives may be incompatible. Consult with pharmacist, if available. When introducing additives, use aseptic technique, mix thoroughly and do not store.

Usual dosage: See package insert.

The overwrap is a moisture barrier.

Use immediately once removed from overwrap.

STORE AT: 20° to 25°C (68° to 77°F) [see USP Controlled Room Temperature]. Protect from freezing.

The container closure is not made with natural rubber latex. Non-PVC, Non-DEHP, Sterile.

BD

LOT

EXP

Becton, Dickinson and Company
1 Becton Drive, Franklin Lakes, NJ 07417 USA
For product inquiry: 1-800-523-0502
Distributed by BD.
Manufactured by Fresenius Kabi.
Made in Germany

403478

0743801.00 US

(01)00317271700073

free flex[®]

500 mL

0.45% Sodium Chloride Injection, USP

For intravenous use. Rx only

2.25 grams per 500 mL (4.5 mg per mL)

Each 100 mL contains: Sodium Chloride 450 mg in water for injection.

Electrolytes per 1,000 mL:

Sodium Chloride 77 mEq 77 mEq

154 mOsmol/LITER (CALC.) pH 5.6 (4.5 to 7.0)

Single-Dose Container. Discard Unused Portion.

Additives may be incompatible. Consult with pharmacist, if available. When introducing additives, use aseptic technique, mix thoroughly and do not store.

Usual dosage: See package insert.

The overwrap is a moisture barrier.

Use immediately once removed from overwrap.

STORE AT: 20° to 25°C (68° to 77°F) [see USP Controlled Room Temperature]. Protect from freezing.

The container closure is not made with natural rubber latex. Non-PVC, Non-DEHP, Sterile.

BD

LOT

EXP

Becton, Dickinson and Company
1 Becton Drive, Franklin Lakes, NJ 07417 USA
For product inquiry: 1-800-523-0502
Distributed by BD.
Manufactured by Fresenius Kabi.
Made in Germany

403477

0743631.00 US

(01)00317271700066

0.45% Sodium Chloride Injection, USP 1000mL

NDC 17271-700-07

SKU: 1727170007

0.45% Sodium Chloride Injection, USP 500mL

NDC 17271-700-06

SKU: 1727170006

free flex[®]

250 mL

0.45% Sodium Chloride Injection, USP

For intravenous use. Rx only

1.125 grams per 250 mL (4.5 mg per mL)

Each 100 mL contains: Sodium Chloride 450 mg in water for injection.

Electrolytes per 1,000 mL:

Sodium Chloride 77 mEq 77 mEq

154 mOsmol/LITER (CALC.) pH 5.6 (4.5 to 7.0)

Single-Dose Container. Discard Unused Portion.

Additives may be incompatible. Consult with pharmacist, if available. When introducing additives, use aseptic technique, mix thoroughly and do not store.

Usual dosage: See package insert.

The overwrap is a moisture barrier.

Use immediately once removed from overwrap.

STORE AT: 20° to 25°C (68° to 77°F) [see USP Controlled Room Temperature]. Protect from freezing.

The container closure is not made with natural rubber latex. Non-PVC, Non-DEHP, Sterile.

BD

LOT

EXP

Becton, Dickinson and Company
1 Becton Drive, Franklin Lakes, NJ 07417 USA
For product inquiry: 1-800-523-0502
Distributed by BD.
Manufactured by Fresenius Kabi.
Made in Norway.

1 234 567 890

403512

FDE 1342 01-69-19-013

(01)00317271700059

free flex[®]

100 mL

0.45% Sodium Chloride Injection, USP

For intravenous use. Rx only

0.45 grams per 100 mL (4.5 mg per mL)

Each 100 mL contains: Sodium Chloride 450 mg in water for injection.

Electrolytes per 1,000 mL:

Sodium Chloride 77 mEq 77 mEq

154 mOsmol/LITER (CALC.) pH 6.0 (4.5 to 7.0)

Single-Dose Container. Discard Unused Portion.

Additives may be incompatible. Consult with pharmacist, if available. When introducing additives, use aseptic technique, mix thoroughly and do not store.

Usual dosage: See package insert.

The overwrap is a moisture barrier.

Use immediately once removed from overwrap.

STORE AT: 20° to 25°C (68° to 77°F) [see USP Controlled Room Temperature]. Protect from freezing.

The container closure is not made with natural rubber latex. Non-PVC, Non-DEHP, Sterile.

BD

LOT

EXP

Becton, Dickinson and Company
1 Becton Drive, Franklin Lakes, NJ 07417 USA
For product inquiry: 1-800-523-0502
Distributed by BD.
Manufactured by Fresenius Kabi.
Made in Norway.

403511

1 234 567 890

FDE 1341 01-69-19-012

(01)00317271700035

0.45% Sodium Chloride
Injection, USP 250mL
NDC 17271-700-05
SKU: 1727170005

freeflex[®]

NDC 17271-700-02

50 mL

0.45% Sodium Chloride Injection, USP

0.225 grams per 50 mL

(4.5 mg per mL)

For intravenous use.

Rx only

Each 100 mL contains: Sodium Chloride

450 mg in water for injection.

Electrolytes per 1,000 mL:

Sodium

77 mEq

Chloride

77 mEq

154 mOsmol/LITER (CALC.)

pH 6.0 (4.5 to 7.0)

Single-Dose Container. Discard Unused Portion.

Additives may be incompatible. Consult with pharmacist, if available. When introducing additives, use aseptic technique, mix thoroughly and do not store.

Usual dosage: See package insert.

The overwrap is a moisture barrier.

Use immediately once removed from overwrap.

STORE AT: 20° to 25°C (68° to 77°F) [see USP Controlled Room Temperature]. Protect from freezing.

The container closure is not made with natural rubber latex.

Non-PVC, Non-DEHP, Sterile.

BD

Becton, Dickinson and Company

1 Becton Drive, Franklin Lakes,

NJ 07417 USA

For product inquiry: 1-800-523-0502

Distributed by BD.

Manufactured by Fresenius Kabi.

Made in Norway.

LOT

EXP

1234567890

403510

FDE 1340 01-69-19-011

(01)00317271700028

0.45% Sodium Chloride
Injection, USP 100mL
NDC 17271-700-03
SKU: 1727170003

The image shows two medical infusion devices. On the left is the BD Alaris™ PCU (Patient Controlled Unit), which has a large color touchscreen displaying a patient interface with buttons for 'NEW PATIENT?', 'STOP/PAUSE/RESTART', and 'Select Yes or No'. Below the screen is a numeric keypad. On the right is the BD Alaris™ Pump, a smaller device with a digital display showing '0.00' and '1.00', and several control buttons. Both devices are white and blue.

0.45% Sodium Chloride
Injection, USP 50mL
NDC 17271-700-02
SKU: 1727170002

BD Alaris™ LVP BUR 60D
LOW SORB 2SS 0.2M BV
SKU: 11613191



**BD Alaris™ LVP 20D LOW
SORB 2SS 0.2M CV**

AS LVP 20D LOW SORB 2SS 0.2M CV
SKU: 11532269



**BD Alaris™ LVP 20D 3SS CV
BV**

AS LVP 20D 3SS CV BV
SKU: 11522558



BD Alaris™ LVP 20D CV

AS LVP 20D 5 SS CV
SKU: 11426965



**BD Alaris™ LVP 20D LOW
SORB**

AS LVP 20D LOW SORB
SKU: 11426964



BD Alaris™ LVP 20D MANI 3WSPK CV

AS LVP 20D MANI 3WSPK CV

SKU: 11419365



BD Alaris™ LVP 10D 3SS CV

SKU: 11404930



Alaris™ Pump Module Set

SKU: 11171447-04



BD Alaris™ LVP 20D 3SS 2CV

AS LVP 20D DEHP SS CV

SKU: 11171447



Alaris™ Pump Module Set

SKU: 10942011-04



BD Alaris™ LVP 20D

AS LVP 20D

SKU: 10942011



BD Alaris™ LVP 20D 4SS 4SS 3G-4WSPK CV

AS LVP 20D 4SS 4SS 3G-4WSPK CV
SKU: 10813621



Gravity administration set

GRAV SET 110" 10DRP V/NV BCV
3NFC-SS ROLL CLMP 2-PIECE LL
SKU: 10802688




BD Alaris™ PCA KIT ASV MICROBORE

PCA KIT ASV MICROBORE
SKU: 10800176



BD Alaris™ PCA KIT ASV MICROBORE CV

PCA KIT ASV MICROBORE CV
SKU: 10800175



Alaris™ SYR PSD ASV 20D SMBORE SS CV

AS SYR PSD ASV 20D SMBORE SS CV
SKU: 10798703



BD Alaris™ SYR PSD AMBER MICROBORE

AS SYR PSD AMBER MICROBORE
SKU: 10798697



Gravity administration set

SKU: 10073124



**BD Alaris™ LVP BLD180M
15D**

AS LVP BLD180M 15D DEHP
SKU: 10062818



Gravity administration set

SKU: 10053857



**BD Alaris™ LVP 10D 3SS 2G-
4WSPK CV**

AS LVP 10D 3SS 2G-4WSPK CV
SKU: 10015896



**BD Alaris™ LVP SS BAG 20D
LOW SORB**

AS LVP SS BAG 20D LOW SORB
SKU: 10015862



**BD Alaris™ LVP SS BAG 20D
LOW SORB 3SS CV**

AS LVP SS BAG 20D LOW SORB 3SS
CV
SKU: 10015861A



**BD Alaris™ LVP BLD180M
15D DEHP**

SKU: 10015414



**BD Alaris™ LVP BUR 20D
LOW SORB SMBORE SS
0.2M**

AS LVP BUR 20D LOW SORB
SMBORE SS 0.2M
SKU: 10015012



**BD Alaris™ SYR PSD LOW
SORB MICROBORE**

AS SYR PSD LOW SORB MICROBORE
SKU: 10014918

**Alaris™ SYR PSD
MICROBORE 0.2M**

AS SYR PSD MICROBORE 0.2M
SKU: 10014916



**Alaris™ SYR PSD
MICROBORE**

AS SYR PSD MICROBORE
SKU: 10014914



**BD Alaris™ LVP 20D LOW
SORB**

SKU: 10013890



**BD Alaris™ LVP SS BAG 20D
2SS TEX CV**

AS LVP SS BAG 20D 2SS TEX CV
SKU: 10013361T



**BD Alaris™ LVP SS BAG 20D
2SS CV**

AS LVP SS BAG 20D 2SS CV
SKU: 10013361



BD Alaris™ LVP 60D DEHP 3SS CV


SKU: 10013186



BD Alaris™ LVP 20D LOW SORB 2SS CV

AS LVP 20D LOW SORB 2SS CV

SKU: 10013072



BD Alaris™ LVP 20D DEHP SS CV

SKU: 10012645



Alaris™ MS III 20D DEHP LOW SORB 2SS CV

SKU: 10011433



BD Alaris™ LVP 20D DEHP 3SS 15M CV

SKU: 10561554

BD Alaris™ LVP SS BAG 20D TEX

SKU: 10321213T



Gravity administration set

SKU: 10010903



BD Alaris™ LVP VSA SMBORE

AS LVP VSA SMBORE

SKU: 10010483



BD Alaris™ LVP 20D LOW SORB SS 0.2M

AS LVP 20D LOW SORB SS 0.2M

SKU: 10010454



BD Alaris™ LVP 20D LOW SORB SS 1.2M

AS LVP 20D LOW SORB SS 1.2M

SKU: 10010453



По вопросам продаж и поддержки обращайтесь:

Алматы (727)345-47-04
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +(727)345-47-04

Беларусь +(375)257-127-884

Узбекистан +998(71)205-18-59

Киргизия +996(312)96-26-47

эл.почта: bog@nt-rt.ru || сайт: <https://becton.nt-rt.ru/>